

Certificate of Quality

VWR Catalogue Number	See List
Description	VWR® Dissecting Forceps, Medium Tip, Serrated
Country of Origin	Manufactured in Pakistan
Date of Issue (yyyy-mm-dd)	2026-06-15

VWR Catalogue Number	Length	Unit
82027-388	122 mm (4.8")	Each
		Pack of 12
		Case of 60
82027-398	114 mm (4½")	Each
		Pack of 12
		Case of 60

Quality System Compliance

This certificate confirms that the product has been sourced from a supplier or manufacturer who has declared their Quality Management System (QMS) complies with the [ISO 9001:2015 & ISO 13485:2016 & F-899](#) standard.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable product specifications.

Product Specifications

Material	AISI 410.
Sterilization	Not applicable.
ATP Assay	Not applicable.
DNase & RNase Free	Not applicable.

BSE/TSE	No use of any raw material produced from, or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
Cytotoxicity	Not applicable.
Non-Pyrogenic Statement	Not applicable.
Latex Statement	The products are free from Latex.
BPA Statement	Not applicable.
DEHP Statement	Not applicable.
RoHS	Not applicable.
REACH Statement	Not applicable.
Storage Conditions	Store in dry place, avoid direct sunlight.
Shelf Life	5 years from date of manufacture.

Disclaimer: VWR'S Quality Management System (QMS) is ISO 9001:2015 certified as a wholesale distributor of finished products, ensuring effective oversight of finished product handling, traceability and distribution processes. VWR is not a manufacturer and does not engage in manufacturing, co- manufacturing, or processing activities such as design, packaging, labelling, testing, or quality release. Consequently, the supplier or manufacturer's QMS certification may differ from VWR's.

VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. This product is intended for research use only. It is not for use in diagnostic procedures, therapeutic applications, or any human or clinical use. This document has been produced electronically and is valid without a signature.

Version History

Sl.no	Version number	Description	Month & Year
1	1.0	Created General CoQ using new layout	March 2021
2	2.0	Updated Quality System Compliance, Shelf life, Disclaimer statement	November 2025
3	2.1	Added new product to the ordering table	June 2026