

Certificate of Conformity

VWR Catalogue Number	See List
Description	Biohazard Bags, Polypropylene, Non-Sterile
Country of Origin	Manufactured in USA
Date of Issue (yyyy-mm-dd)	2024-08-23

Catalogue Number	Capacity [L]	Dimension [mm]	Thickness [µm]
129-0030	12	610x305	75
129-0031	60	914x610	75
129-0032	200	1220x965	75
129-1560	45	762x610	76
129-1561	60	914x610	100
129-1562	100	1016x787	100

Quality System Compliance

Products are manufactured under the **ISO 9001:2015** standard. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

QC Testing

Representative products samples are tested in accordance with internationally recognized ASTM Test methods for Dart Impact Strength (ASTM D1709) and tear resistance (ASTM D1922).

Product Specifications

Material	Polypropylene. Compliant with Federal DOT Regulation 49 CFR 173.197.
Autoclavable	The products allows use for autoclave sterilizer operating temperatures of 135 °C.
Sterilization	Not Applicable.
ATP Assay	Not Applicable.

DNase & RNase Free	Not Applicable.
BSE/TSE	Resins and films used to manufacture this product may contain byproducts of bovine sourced raw material in their manufacture. These animal byproduct additives are processed under highly refined and controlled industrial chemical processes and at very high temperatures prior to purification and subsequent use as a polymer additive. The animal derived byproducts are not added during the process of manufacturing but originate either in the manufacture of the resins or in the additives used in the film formulations.
Cytotoxicity	Not Applicable.
Non-Pyrogenic Statement	Not Applicable.
Latex Statement	This product is not made from natural rubber latex.
BPA Statement	Not Applicable.
DEHP Statement	Please see RoHS information below.
RoHS	No substances (Lead, Cadmium, Mercury, Hexavalent Chromium (Cr6+), Poly Brominated Biphenyls (PBB), Poly Brominated Diphenyl ethers (PBDE), Bis(2-Ethyhexyl) Phthalate (DEHP), Benzyl Butyl Phthalate (BBP), Dibutyl Phthalate (DBP), Diisobutyl Phthalate (DIBP)) are used in manufacturing the raw materials and final product. No routinely analyse is performed.
REACH Statement	No substances as listed on >Substance of Very High Concern< candidate list published 2022-01-17 are used in manufacturing of the raw materials used to produce this product. No routinely analysis is performed.
Storage Conditions	Room temperature.
Shelf Life	No Shelf life.

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
This document has been produced electronically and is valid without a signature.

Annex

TYPICAL PRODUCT CHARACTERISTICS:

Properties		2 mil (50.8μ)	3 mil (76μ)	4 mil (102μ)	Units	Test Method
Elmendorf Tear Resistance*	MD ↕ TD ↔	600 900	600 900	600 900	grams grams	ASTM D1922
Dart Impact Strength		>1400	>1000	>800	grams	ASTM D1709
Coefficient of Friction		<0.20	<0.20	<0.20		ASTM D1894
Yield		15300	10200	7650	sq. in./lb.	Internal