

Certificate of Quality

VWR Catalogue Number	89132-060
Description	VWR® Square PETG Sterile Media Bottle with stranded caps, Capacity 250 mL (8.5 oz.), Cap size 38-430
Lot Number	111013Q
Date of Manufacture (yyyy-mm-dd)	2022-10-13
Date of Sterilization (yyyy-mm-dd)	2022-12-13
Date of Expiration (yyyy-mm-dd)	4 years post sterilization
Country of Origin	Manufactured in Taiwan
Date of Issue (yyyy-mm-dd)	2023-09-19

Quality System Compliance

Products are manufactured under the ISO 9001:2015 standard. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements.

Product Specifications

This product was manufactured in a class 10,000 cleanroom using FDA compliant SK SKYGREEN S2008 resin. This product is deemed particle free.

Material

Polyethylene terephthalate glycol (PETG)
Height: 141 mm; Neck O.D.: 34.5 mm; Outer diameter: 61 mm

All component materials have been tested and met the requirements for United States Pharmacopoeia (USP) Class VI Biological Test for Plastics, current edition.

The plastics meet the requirements of the United States Food and Drug Administration (FDA) for food and beverage contact in 21 CFR 177.1640.



中國生化科技股份有限公司

CHINA BIOTECH CORPORATION

誠信正直 · 顧客感動 · 社會公道

TEL: 886-4-23597515 FAX: 886-4-23597080
台中市工業區33路10號
10, 33rd Road, Taichung Industrial Park,
Taichung, Taiwan R.O.C 407

照射證明書
CERTIFICATE OF IRRADIATION

DATE: 2022/12/20

行政院原能會核准設立照射廠執照証號 IRRADIATION PLANT NO: 物字第 1100223 號

客戶名稱 CUSTOMER NAME: TRIFOREST ENTERPRISES, INC.

照射日期 IRRADIATION RUN DATE: 2022/11/26

照射批號 IRRADIATION RUN NUMBER: SAJ22086-J

客戶產品已照射 MATERIALS PROCESSED:

箱數 CTNS	內容 DESCRIPTION	客戶產品批號 LOT NO
83	250 ml PETG SQUARE MEDIA Bottle, Sterile (BGC0250SB-VWR)	111013Q
40	250 ml PETG SQUARE MEDIA Bottle, Sterile (BGC0250S-TF)	111018Q
23	250 ml PETG SQUARE MEDIA Bottle, Sterile (BGC0250SB-VWR)	111108Q
8	250 ml PETG SQUARE MEDIA Bottle, Sterile (BGC0250SB-VWR)	111102Q
總數 TOTAL	154 箱數 CTNS	

中國生化科技股份有限公司證明上述產品經本公司劑量偵測系統判讀，吸收劑量如下：
China Biotech Corporation certifies that the material listed above (has described by its manufacturer)
received the following doses within the precision limits of the dosimetry system employed

最低劑量 23.3 kGy; 最高劑量 27.7 kGy
MINIMUM DOSAGE MAXIMUM DOSAGE

使用放射性同位素 ISOTOPE UTILIZED: 鈷 60 COBALT-60

客戶劑量要求 DOSE REQUIREMENT: 最低劑量 MIN 18.0 kGy; 最高劑量 MAX 30.0 kGy

確認者:

CERTIFIED BY

品質部主管

QUALITY ASSURANCE