

BAXTER CLEVELAND CERTIFICATE OF ANALYSIS

911 North Davis
Cleveland, MS 38732

(2F7121 / 2F7122 / 2F7123 / 2F7124 / JF7123 / JF7124)

CODE 2F7123	BATCH M162654	MFG. DATE 3-10-23	STERILIZATION DATE 3-10-23
MANUFACTURING FACILITY CLEVELAND, MS		EXPIRY DATE Mar 26	

DESCRIPTION

0.9% Sodium Chloride, USP

TEST	LIMITS	SAMPLE RESULTS	DATE PASSED
BACTERIAL ENDOTOXIN (PYROGEN)	CONTROL LIMIT (Eu/mL)		
13-01-10-313	<0.5 EU/ml = Pass	EU/ml 20.01	3-13-23
STERILITY PARAMETRIC RELEASE	PASS	PASSED	3-10-23
PARTICLE ANALYSIS 20-03-07-006	NMT 25 part/mL $\geq 10 \mu\text{m}$	part/mL 1	3-14-23
	NMT 3 part/mL $\geq 25 \mu\text{m}$	part/mL 0	

PARAMETRIC RELEASE IN ACCORDANCE WITH U.S. FDA REQUIREMENTS

CHEMICAL RESULTS 3-15-19-156 ISSUE DATE **11-13-18**

REFERENCE	TEST	LIMITS	RESULT
11-21-16-003	pH @ 25°C	4.5-7.0	5.6
11-07-19-001 or 11-16-03-005 or 11-16-03-016	NaCl Assay	8.55-9.45 g/l	8.89 g/l
11-20-19-006	Sodium ID - USP	Positive	Positive
			NA
			E. 3-20-23

THIS PRODUCT WAS MANUFACTURED ACCORDING TO CURRENT GMP, FULLFILLS THE REQUIRMENTS OF THE MASTER
PRODUCTION DOCUMENT AND TESTED TO USP COMPLIANCE.

DATE BATCH RELEASED **3-20-23**PREPARED BY/DATE **E. Solis 3-20-23**QUALITY MANAGEMENT SUPERVISOR **C. Hogan 3-22-23**

MDM

1. Code, Batch, Mfg. Date, Sterilization Date, and Expiry Date - Enter the information from Production Report of batch required.

2. Test:

Bacterial Endotoxin – Enter the bacterial endotoxin test results. Enter the date passed (date the bacterial endotoxin testing was complete) from the bacterial endotoxin test report for the batch.

Sterility – Enter the date product was sterilized.

Particle Analysis – Enter the sample results and the date passed from the particle analysis record for the batch.

3. Chemical Results - Enter issue date of the applicable 3-15 specification. All results should be verified against the limits in the applicable 3-15 specification. All test results are retrieved from the EOB chemical testing. Fill in test results.

Note: Multiple test results may be averaged by Release Clerk and the average entered on this certification.
(Example - multitank and stability batches).

Note: NaCl Assay- Circle the 11 spec used for testing.

4. Date Batch Released - Fill in date batch was released.
5. Prepared By/Date - Dated signature of person entering information.

Note: Signature of Release Clerk indicates the results have been verified against the results listed on the chemical final analysis and the applicable 3-15 chemical testing specification, the results listed on the Particle Analysis Report and results listed on the Bacterial Endotoxin report.

6. N/A sections not used.

7. Quality Management Supervisor - Signature of Quality Management Supervisor or delegate verifying the Batch #, Code, Manufacturer Date, Expiration Date, and Batch Release Date documented matches the information documented on the Production Record. Individual will and verify that all test results listed are within specification and/or SOP limits.

DESCRIPTION OF CHANGE		
Name change. Update title.		
REASON FOR CHANGE		
Update approvals on file. Update title of form.		
CHANGE HISTORY		
CHANGE #	ISSUE DATE	DESCRIPTION OF CHANGE
21-1779	01/11/22	Name change.
20-2021	07/16/20	Add JF7123 and JF7124 to form CLF-03-148.
20-1012	04/16/20	Update MDM to include specific information for the approver to verify.
APPROVALS ON FILE		
ROLE		NAME
SME/Release Clerk		Beverly Robinson
Approver		Curtissha Hogan
Concurrer		Ellen Soliz