



Product Specification Sheet

VWR MICROSCOPE SLIDES, GROUND EDGES 45°

VWR International certifies that the specified products meet all VWR specifications and quality criteria for release. The products comply in every respect with the relevant particular specifications, drawings, relevant standards and regulations in force.

VWR International certifie que les produits mentionnés respectent toutes les spécifications et critères de qualité VWR. Les produits satisfont aux spécifications particulières, aux plans, aux normes et règlements en vigueur s'y rapportant.

VWR International bestätigt hiermit, daß die angegebenen Produkte alle VWR Spezifikationen und Qualitätskriterien zur Freigabe erfüllen. Die Produkte entsprechen in jeder Hinsicht den maßgeblichen Spezifikationen, den Zeichnungen sowie den relevanten Normen und gültigen Vorschriften.

Cat. No.	Colour	Pk
631-1560	White	50
631-1561	Pink	50
631-1562	Blue	50
631-1563	Yellow	50
631-1564	Green	50
631-1565	Orange	50

Material: Clear white glass

W×D×H: 76×26×1 mm

Other characteristics:

- Frosted ends in six different colours
- Writing surface is resistant to all common dyes and solvents

Disclaimer: VWR wants to state that these specifications can only be regarded as general information and will not discharge the user from his obligation to ensure the product is suitable for the intended use - If there is any doubt ask us.

Date of issue: **June 2013**