




100-22027730


Description	MEL Number	Sample vol μl	Quantity(PCS) per pack	Pack per case	Cover Type	Path length mm	Material	A	B	C	D	E	F	G	Sides Color
Absorption Quartz cuvettes with flat lid and round bottom	VWRI634-1196	350	1	8	Flat PTFE LID	1	Quartz Glass	12.5±0.1	10±0.2	3.5±0.1	1±0.03	45±0.5	43.5±0.5	1.25	2 Frosted Sides
	VWRI634-1197	700	1	8	Flat PTFE LID	2	Quartz Glass	12.5±0.1	10±0.2	4.5±0.1	2±0.03	45±0.5	43.5±0.5	1.25	2 Frosted Sides
	VWRI634-1198	1750	1	8	Flat PTFE LID	5	Quartz Glass	12.5±0.1	10±0.2	7.5±0.1	5±0.03	45±0.5	43.5±0.5	1.25	2 Frosted Sides
	VWRI634-1199	3500	1	8	Flat PTFE LID	10	Quartz Glass	12.5±0.1	10±0.2	12.5±0.1	10±0.03	45±0.5	43.5±0.5	1.25	2 Frosted Sides
	VWRI634-1200	3500	2	8	Flat PTFE LID	10	Quartz Glass	12.5±0.1	10±0.2	12.5±0.1	10±0.03	45±0.5	43.5±0.5	1.25	2 Frosted Sides
	VWRI634-1201	7000	1	8	Flat PTFE LID	20	Quartz Glass	12.5±0.1	10±0.2	22.5±0.1	20±0.03	45±0.5	43.5±0.5	1.25	2 Frosted Sides
	VWRI634-1202	14000	1	8	Flat PTFE LID	40	Quartz Glass	12.5±0.1	9.5±0.2	42.5±0.1	40±0.05	45±0.5	43.5±0.5	1.25	2 Frosted Sides
	VWRI634-1203	17500	1	8	Flat PTFE LID	50	Quartz Glass	12.5±0.1	9.5±0.2	52.5±0.1	50±0.05	45±0.5	43.5±0.5	1.25	2 Frosted Sides
	VWRI634-1204	35000	1	4	Flat PTFE LID	100	Quartz Glass	12.5±0.1	9.5±0.2	102.5±0.1	100±0.05	45±0.5	43.5±0.5	1.25	2 Frosted Sides



SECTION X - X



SECTION Y - Y



FROSTED



Notes:

- Product data: see product data table
- Product standard: Q/320282AAE001-2020
- Matched cells exhibit a transmission tolerance within 0.5% at 200nm
- Transmission: > 80% (at 200nm)
- Label: comply with VWR label guideline-issued separately
- Carton box: comply with VWR packing requirements-issue separately
- Major defects: AQL 1.0 per General Inspection Plan Level II, single sampling (ANSI/ASQC Z1.4-1993)
 - Box is damaged,punctured and not usable
 - Cuvettes appear dirty or damaged
 - Label or printing on label missing, illegible or incorrect
 - Transmission doesn't meet spec.
 - Dimension is out of tolerance
- First Article Report required on new product.
- Supplier is required to inform VWR and obtain approval before any changes to materials or source, process change or manufacturing location change.
- Must be manufactured at supplier's premises condition and not from 3rd party



MARK	ECO NO.		ZONE	EDITED	DATE	APPD	DATE	REV
APPROVALS	BY	DATE	THIS DOCUMENT CONTAINS CONFIDENTIAL AND PROPRIETARY INFORMATION AND MAY NOT BE COPIED OR REPRODUCED IN WHOLE OR IN PART NOR PUT ON CIRCULATIONS WITHOUT THE EXPRESS WRITTEN PERMISSION OF VWR INTERNATIONAL, INC.					
DRAWN		8-Apr-2022						
DESIGN			UNLESS OTHERWISE SPECIFIED TOLERANCES AS FOLLOWS DIMENSIONS X ± 0.5 (XX ± 0.25) - DIMENSIONS ARE IN MM UNLESS XX ± 0.25 (XXX ± 0.1) - SURFACES FINISH TO BE XXX ± 0.2 (XXXX ± 0.05) - DENOTES CHANGE TO BE ANGLES X ± 1° (XX ± 0.5°) - WORKMANSHIP SHALL BE THE XX ± 0.5° BEST COMMERCIAL QUALITY					
MFG								
CHECKED			FINISH: N/A					
APPROVED								
			VWR P/N		REV. A0	SIZE: A3		
			REFERENCE DRAWING NO.		SCALE: N/A			
			MATERIAL:		WEIGHT:	SHEET 1 OF 1		

200-22027730

Description	MEL Number	Sample vol μl	Quantity(PCS) per pack	Pack per case	Cover Type	Path length mm	Material	A	B	C	D	E	F	G	Sides Color
Absorption optical cuvettes with flat lid and round bottom	VWRI634-1205	3500	1	8	Flat PTFE LID	10	Optical Glass	12.5±0.1	10±0.2	12.5±0.1	10±0.03	45±0.5	43.5±0.5	1.25	2 Frosted Sides
	VWRI634-1206	3500	2	8	Flat PTFE LID	10	Optical Glass	12.5±0.1	10±0.2	12.5±0.1	10±0.03	45±0.5	43.5±0.5	1.25	2 Frosted Sides
	VWRI634-1207	7000	1	8	Flat PTFE LID	20	Optical Glass	12.5±0.1	10±0.2	22.5±0.1	20±0.03	45±0.5	43.5±0.5	1.25	2 Frosted Sides
	VWRI634-1208	14000	1	8	Flat PTFE LID	40	Optical Glass	12.5±0.1	9.5±0.2	42.5±0.1	40±0.05	45±0.5	43.5±0.5	1.25	2 Frosted Sides
	VWRI634-1209	17500	1	8	Flat PTFE LID	50	Optical Glass	12.5±0.1	9.5±0.2	52.5±0.1	50±0.05	45±0.5	43.5±0.5	1.25	2 Frosted Sides
	VWRI634-1210	35000	1	4	Flat PTFE LID	100	Optical Glass	12.5±0.1	9.5±0.2	102.5±0.1	100±0.05	45±0.5	43.5±0.5	1.25	2 Frosted Sides



SECTION X - X



SECTION Y - Y


Notes:

- 1.Product data: see product data table
- 2.Product standard: Q/320282AAE001-2020
- 3.Matched cells exhibit a transmission tolerance within 0.5% at 360nm
- 4.Transmission: > 80% (at 360nm)
- 5.Label: comply with VWR label guideline-issued separately
- 6.Carton box: comply with VWR packing requirements-issue separately
- 7.Major defects: AQL 1.0 per General Inspection Plan Level II, single sampling (ANSI/ASQC Z1.4-1993)
 - a. Box is damaged,punctured and not usable
 - b. Cuvettes appear dirty or damaged
 - c. Label or printing on label missing, illegible or incorrect
 - d. Transmission doesn't meet spec.
 - e. Dimension is out of tolerance
8. First Article Report required on new product.
9. Supplier is required to inform VWR and obtain approval before any changes to materials or source, process change or manufacturing location change.
10. Must be manufactured at supplier's premises condition and not from 3rd party


MARK	ECO NO.		ZONE	EDITED	DATE	APPD	DATE	REV
APPROVALS	BY	DATE	THIS DOCUMENT CONTAINS CONFIDENTIAL AND PROPRIETARY INFORMATION AND MAY NOT BE COPIED OR REPRODUCED IN WHOLE OR IN PART NOR PUT ON CIRCULATIONS WITHOUT THE EXPRESS WRITTEN PERMISSION OF VWR INTERNATIONAL, INC. UNLESS OTHERWISE SPECIFIED TOLERANCES AS FOLLOWS DIMENSIONS X ± 0.5 (XXX ± 0.02) XX ± 0.25 (XXX ± 0.01) XX ± 0.2 (XXX ± 0.005) ANGLES X ± 1° XX ± 5° - DIMENSIONS ARE IN MM UNLESS - SURFACES FINISH TO BE - DENOTES CHANGE - WORKMANSHIP SHALL BE THE BEST COMMERCIAL QUALITY FINISH: N/A					
DRAWN		8-Apr-2022						
DESIGN								
MFG								
CHECKED			DESCRIPTION: Absorption optical cuvettes with flat lid and round bottom VWR P/N: CELL2022-002 REV. A0 SIZE: A3 REFERENCE: DRAWING NO. SCALE: N/A					
APPROVED								
			MATERIAL:		WEIGHT:		SHEET 1 OF 1	

CELL2022-003

Description	MEL Number	Sample vol μl	Quantity(PCS) per pack	Pack per case	Cover Type	Path length mm	Material	A	B	C	D	E	F	G	H	I	Sides Color
Absorption quartz cuvettes with Teflon stopper and flat bottom	VWRI634-1211	350	1	8	PTFE stopper	1	Quartz Glass	12.5±0.1	10±0.2	3.5±0.1	1±0.03	52±0.5	50.5±0.5	1.25	8	Ø3±0.5	2 Frosted Sides
	VWRI634-1212	700	1	8	PTFE stopper	2	Quartz Glass	12.5±0.1	10±0.2	4.5±0.1	2±0.03	52±0.5	50.5±0.5	1.25	8	Ø3±0.5	2 Frosted Sides



SECTION X - X



SECTION Y - Y



Notes:

- Product data: see product data table
- Product standard: Q/320282AAE001-2020
- Transmission: > 80% (at 200nm)
- Label: comply with VWR label guideline-issued separately
- Carton box: comply with VWR packing requirements-issue separately
- Major defects: AQL 1.0 per General Inspection Plan Level II, single sampling (ANSI/ASQC Z1.4-1993)
 - Box is damaged,punctured and not usable
 - Cuvettes appear dirty or damaged
 - Label or printing on label missing, illegible or incorrect
 - Transmission doesn't meet spec.
 - Dimension is out of tolerance
- First Article Report required on new product.
- Supplier is required to inform VWR and obtain approval before any changes to materials or source, process change or manufacturing location change.
- Must be manufactured at supplier's premises condition and not from 3rd party

MARK	ECO NO.		ZONE	EDITED	DATE	APPD	DATE	REV
APPROVALS	BY	DATE	THIS DOCUMENT CONTAINS CONFIDENTIAL AND PROPRIETARY INFORMATION AND MAY NOT BE COPIED OR REPRODUCED IN WHOLE OR IN PART NOR PUT ON CIRCULATIONS WITHOUT THE EXPRESS WRITTEN PERMISSION OF VWR INTERNATIONAL, INC.					
DRAWN		8-Apr-2022						
DESIGN			UNLESS OTHERWISE SPECIFIED TOLERANCES AS FOLLOWS DIMENSIONS X ± 0.5 DXX ± 0.02 XX ± 0.25 DXXX ± 0.01 XXX ± 0.2 DXXXX ± 0.005 ANGLES X ± 1° XX ± 5°					
MFG								
CHECKED			- DIMENSIONS ARE IN MM UNLESS - SURFACES FINISH TO BE $\sqrt{R}$ - DIMENSIONS CHANGE Δ - WORKMANSHIP SHALL BE THE BEST COMMERCIAL QUALITY					
APPROVED								
FINISH: N/A			DESCRIPTION: Absorption quartz cuvettes with Teflon stopper and flat bottom					
			VWR P/N CELL2022-003		REV. A0		SIZE: A3	
			REFERENCE DRAWING NO.				SCALE: N/A	
			MATERIAL:		WEIGHT:		SHEET 1 OF 1	

CELL2022-004

Description	MEL Number	Sample vol μl	Quantity(PCS) per pack	Pack per case	Cover Type	Path length mm	Material	A	B	C	D	E	F	G	H	I	Sides Color
Absorption quartz cuvettes with Teflon stopper and round bottom	VWRI634-1213	1750	1	8	PTFE stopper	5	Quartz Glass	12.5±0.1	10±0.2	7.5±0.1	5±0.03	45±0.5	43.5±0.5	1.25	5	Ø4.5±0.5	2 Frosted Sides
	VWRI634-1214	3500	1	8	PTFE stopper	10	Quartz Glass	12.5±0.1	10±0.2	12.5±0.1	10±0.03	45±0.5	43.5±0.5	1.25	5	Ø9.5±0.5	2 Frosted Sides
	VWRI634-1215	7000	1	8	PTFE stopper	20	Quartz Glass	12.5±0.1	10±0.2	22.5±0.1	20±0.03	45±0.5	43.5±0.5	1.25	5	Ø9.5±0.5	2 Frosted Sides





Notes:

- Product data: see product data table
- Product standard: Q/320282AAE001-2020
- Transmission: > 80% (at 200nm)
- Label: comply with VWR label guideline-issued separately
- Carton box: comply with VWR packing requirements-issue separately
- Major defects: AQL 1.0 per General Inspection Plan Level II, single sampling (ANSI/ASQC Z1.4-1993)
 - Box is damaged,punctured and not usable
 - Cuvettes appear dirty or damaged
 - Label or printing on label missing, illegible or incorrect
 - Transmission doesn't meet spec.
 - Dimension is out of tolerance
- First Article Report required on new product.
- Supplier is required to inform VWR and obtain approval before any changes to materials or source, process change or manufacturing location change.
- Must be manufactured at supplier's premises condition and not from 3rd party



MARK	ECO NO.		ZONE	EDITED	DATE	APPD	DATE	REV
APPROVALS	BY	DATE	THIS DOCUMENT CONTAINS CONFIDENTIAL AND PROPRIETARY INFORMATION AND MAY NOT BE COPIED OR REPRODUCED IN WHOLE OR IN PART NOR PUT ON CIRCULATIONS WITHOUT THE EXPRESS WRITTEN PERMISSION OF VWR INTERNATIONAL, INC.					
DRAWN		8-Apr-2022						
DESIGN			UNLESS OTHERWISE SPECIFIED TOLERANCES AS FOLLOWS DIMENSIONS X ± 0.5 DXX ± 0.02 XX ± 0.25 DXXX ± .01 XXX ± 0.12 DXXXX ± .005 ANGLES X ± 1° XX ± .5°					
MFG								
CHECKED			- DIMENSIONS ARE IN MM UNLESS - SURFACES FINISHES TO BE - DENOTES CHANGE - WORKMANSHIP SHALL BE THE BEST COMMERCIAL QUALITY					
APPROVED								
FINISH: N/A			DESCRIPTION: Absorption quartz cuvettes with Teflon stopper and round bottom					
			VWR P/N: CELL2022-004 REV: A0 SIZE: A3					
			REFERENCE: DRAWING NO. SCALE: N/A					
			MATERIAL: WEIGHT: SHEET 1 OF 1					

CELL2022-005

Description	MEL Number	Sample vol μl	Quantity(PCS) per pack	Pack per case	Cover Type	Path length mm	Material	A	B	C	D	E	F	G	Transmission	Sides Color
Absorption cuvettes with flat lid and flat bottom	VWRI634-1216	3500	1	8	Flat PTFE LID	10	Optical Glass	12.5±0.1	10±0.2	12.5±0.1	10±0.03	45±0.5	43.5±0.5	1.25	360nm> 80%	2 Frosted Sides
	VWRI634-1217	7000	1	8	Flat PTFE LID	20	Optical Glass	12.5±0.1	10±0.2	22.5±0.1	20±0.03	45±0.5	43.5±0.5	1.25	360nm> 80%	2 Frosted Sides
	VWRI634-1218	17500	1	8	Flat PTFE LID	50	Optical Glass	12.5±0.1	9.5±0.2	52.5±0.1	50±0.05	45±0.5	43.5±0.5	1.25	360nm> 80%	2 Frosted Sides
	VWRI634-1219	3500	1	8	Flat PTFE LID	10	Quartz Glass	12.5±0.1	10±0.2	12.5±0.1	10±0.03	45±0.5	43.5±0.5	1.25	200nm> 80%	2 Frosted Sides



SECTION X - X



SECTION Y - Y



Notes:

- Product data: see product data table
- Product standard: Q/320282AAE001-2020
- Label: comply with VWR label guideline-issued separately
- Carton box: comply with VWR packing requirements-issue separately
- Major defects: AQL 1.0 per General Inspection Plan Level II, single sampling (ANSI/ASQC Z1.4-1993)
 - Box is damaged,punctured and not usable
 - Cuvettes appear dirty or damaged
 - Label or printing on label missing, illegible or incorrect
 - Transmission doesn't meet spec.
 - Dimension is out of tolerance
- First Article Report required on new product.
- Supplier is required to inform VWR and obtain approval before any changes to materials or source, process change or manufacturing location change.
- Must be manufactured at supplier's premises condition and not from 3rd party

MARK		ECO NO.		ZONE		EDITED	DATE	APPD	DATE	REV	
APPROVALS	BY	DATE		THIS DOCUMENT CONTAINS CONFIDENTIAL AND PROPRIETARY) INFORMATION AND MAY NOT BE COPIED OR REPRODUCED IN WHOLE OR IN PART NOR PUT ON CIRCULATIONS WITHOUT THE EXPRESS WRITTEN PERMISSION OF VWR INTERNATIONAL, INC UNLESS OTHERWISE SPECIFIED TOLERANCES AS FOLLOWS DIMENSIONS X ± 0.5 DXX ± 0.02 - DIMENSIONS ARE IN MM UNLESS XX ± 0.25 DXXX ± 0.01 - SURFACES FINISHES TO BE XXX ± 0.12 DXXXX ± 0.005 - DENOTES CHANGE Δ ANGLES X ± 1° - WORKMANSHIP SHALL BE THE XX ± 5° BEST COMMERCIAL QUALITY FINISH: N/A		PROJ: 		DESCRIPTION: Absorption cuvettes with flat lid and flat bottom			
DRAWN	8-Apr-2022		VWR P/N CELL2022-005			REV. A0	SIZE: A3				
DESIGN			REFERENCE DRAWING NO.			SCALE: N/A					
MFG			MATERIAL:			WEIGHT:	SHEET 1 OF 1				
CHECKED											
APPROVED											

900-22027730

Description	MEL Number	Sample vol μl	Quantity(PCS) per pack	Pack per case	Cover Type	Path length mm	Material	A	B	C	D	E	F	G	Sides Color
Micro cuvettes with flat lid and flat bottom	VWRI634-1220	1000	1	8	Flat PTFE LID	10	Quartz Glass	12.5±0.1	3±0.2	12.5±0.1	10±0.03	45±0.5	30°	1.25	2 Frosted Sides
	VWRI634-1221	1400	1	8	Flat PTFE LID	10	Quartz Glass	12.5±0.1	4±0.2	12.5±0.1	10±0.03	45±0.5	30°	1.25	2 Frosted Sides
	VWRI634-1222	7000	1	8	Flat PTFE LID	50	Quartz Glass	12.5±0.1	4±0.2	52.5±0.1	50±0.05	45±0.5	30°	1.25	2 Frosted Sides






Notes:

- 1.Product data: see product data table
- 2.Product standard: Q/320282AAE001-2020
- 3.Transmission: > 80% (at 200nm)
- 4.Label: comply with VWR label guideline-issued separately
- 5.Carton box: comply with VWR packing requirements-issue separately
- 6.Major defects: AQL 1.0 per General Inspection Plan Level II, single sampling (ANSI/ASQC Z1.4-1993)
 - a. Box is damaged,punctured and not usable
 - b. Cuvettes appear dirty or damaged
 - c. Label or printing on label missing, illegible or incorrect
 - d. Transmission doesn't meet spec.
 - e. Dimension is out of tolerance
7. First Article Report required on new product.
8. Supplier is required to inform VWR and obtain approval before any changes to materials or source, process change or manufacturing location change.
9. Must be manufactured at supplier's premises condition and not from 3rd party


MARK	ECO NO.		ZONE		EDITED	DATE	APPD	DATE	REV
APPROVALS	BY	DATE	THIS DOCUMENT CONTAINS CONFIDENTIAL AND PROPRIETY) INFORMATION AND MAY NOT BE COPIED OR REPRODUCED IN WHOLE OR IN PART NOR PUT ON CIRCULATIONS WITHOUT THE EXPRESS WRITTEN PERMISSION OF VWR INTERNATIONAL, INC.		FIG. 1 				
DRAWN		8-Apr-2022							
DESIGN									
MFG									
CHECKED									
APPROVED			UNLESS OTHERWISE SPECIFIED TOLERANCES AS FOLLOWS DIMENSIONS X ± 0.5 DXX ± 0.02 - DIMENSIONS ARE IN MM UNLESS XX ± 0.25 DXXX ± 0.01 - SURFACES FINISH TO BE $\sqrt{\text{ }}$ XXX ± 0.2 DXXXX ± 0.05 - DENOTES CHANGE Δ ANGLES X ± 1° - WORKMANSHIP SHALL BE THE XX ± 5° BEST COMMERCIAL QUALITY		DESCRIPTION: Micro cuvettes with flat lid and flat bottom				
FINISH: N/A			VWR P/N		CELL2022-006		REV. A0	SIZE: A3	
			REFERENCE DRAWING NO.		SCALE: N/A				
			MATERIAL:				WEIGHT:	SHEET 1 OF 1	

CELL2022-007

Description	MEL Number	Sample vol μl	Quantity(PCS) per pack	Pack per case	Cover Type	Path length mm	Material	A	B	C	D	E	F	G	Sides Color
Fluorometer cuvettes with flat lid and flat bottom	VWRI634-1223	3500	1	8	Flat PTFE LID	10	Quartz Glass	12.5±0.1	10±0.03	12.5±0.1	10±0.03	45±0.5	43.5±0.5	1.25	4 Polished Sides
	VWRI634-1224	7000	1	8	Flat PTFE LID	20	Quartz Glass	12.5±0.1	10±0.03	22.5±0.1	20±0.03	45±0.5	43.5±0.5	1.25	4 Polished Sides



SECTION X - X



SECTION Y - Y



Notes:

- Product data: see product data table
- Product standard: Q/320282AAE001-2020
- Transmission: > 80% (at 200nm)
- Label: comply with VWR label guideline-issued separately
- Carton box: comply with VWR packing requirements-issue separately
- Major defects: AQL 1.0 per General Inspection Plan Level II, single sampling (ANSI/ASQC Z1.4-1993)
 - Box is damaged,punctured and not usable
 - Cuvettes appear dirty or damaged
 - Label or printing on label missing, illegible or incorrect
 - Transmission doesn't meet spec.
 - Dimension is out of tolerance
- First Article Report required on new product.
- Supplier is required to inform VWR and obtain approval before any changes to materials or source, process change or manufacturing location change.
- Must be manufactured at supplier's premises condition and not from 3rd party

MARK		ECO NO.		ZONE		EDITED	DATE	APPD	DATE	REV	
APPROVALS	BY	DATE	THIS DOCUMENT CONTAINS CONFIDENTIAL AND PROPRIETY) INFORMATION AND MAY NOT BE COPIED OR REPRODUCED IN WHOLE OR IN PART NOR PUT ON CIRCULATIONS WITHOUT THE EXPRESS WRITTEN PERMISSION OF VWR INTERNATIONAL, INC. UNLESS OTHERWISE SPECIFIED TOLERANCES AS FOLLOWS DECIMALS X ± 0.5 XXX ± 0.02 - DIMENSIONS ARE IN MM UNLESS XXX ± 0.25 XXXX ± 0.01 - SURFACES FINISHES TO BE  XXXX ± 0.2 XXXX ± 0.05 - DENOTES CHANGE Δ ANGLES X ± 1° - WORKMANSHIP SHALL BE THE XX ± 5° BEST COMMERCIAL QUALITY								
DRAWN		8-Apr-2022									
DESIGN											
MFG											
CHECKED											
APPROVED			FINISH: N/A			DESCRIPTION: Fluorometer cuvettes with flat lid and flat bottom					
						VWR P/N	CELL2022-007		REV. A0	SIZE: A3	
						REFERENCE DRAWING NO.				SCALE: N/A	
						MATERIAL:				WEIGHT:	SHEET 1 OF 1

800-22027730

Description	MEL Number	Sample vol μl	Quantity(PCS) per pack	Pack per case	Cover Type	Path length mm	Material	A	B	C	D	E	F	G	H	I	Sides Color
Fluorometer cuvettes with Teflon stopper and flat bottom	VWRI634-1225	3500	1	8	PTFE stopper	10	Quartz Glass	12.5±0.1	10±0.03	12.5±0.1	10±0.03	45±0.5	43.5±0.5	1.25	5	Ø9.5±0.5	4 Polished Sides






Notes:


- Product data: see product data table
- Product standard: Q/320282AAE001-2020
- Transmission: > 80% (at 200nm)
- Label: comply with VWR label guideline-issued separately
- Carton box: comply with VWR packing requirements-issue separately
- Major defects: AQL 1.0 per General Inspection Plan Level II, single sampling (ANSI/ASQC Z1.4-1993)
 - Box is damaged,punctured and not usable
 - Cuvettes appear dirty or damaged
 - Label or printing on label missing, illegible or incorrect
 - Transmission doesn't meet spec.
 - Dimension is out of tolerance
- First Article Report required on new product.
- Supplier is required to inform VWR and obtain approval before any changes to materials or source, process change or manufacturing location change.
- Must be manufactured at supplier's premises condition and not from 3rd party

MARK	ECO NO.				ZONE	EDITED	DATE	APPD	DATE	REV
APPROVALS	BY	DATE	THIS DOCUMENT CONTAINS CONFIDENTIAL AND PROPRIETY) INFORMATION AND MAY NOT BE COPIED OR REPRODUCED IN WHOLE OR IN PART NOR PUT ON CIRCULATIONS WITHOUT THE EXPRESS WRITTEN PERMISSION OF VWR INTERNATIONAL, INC			FIELD 				
DRAWN		8-Apr-2022					DESCRIPTION: Fluorometer cuvettes with Teflon stopper and flat bottom			
DESIGN			UNLESS OTHERWISE SPECIFIED TOLERANCES AS FOLLOWS DECIMALS X ± 0.5 (XX ± 0.02) XX ± 0.25 (XXX ± 0.01) XXX ± 0.12 (XXXX ± 0.005) ANGLES X ± 1° XX ± 3°			VWR P/N	REV. A0	SIZE: A3		
MFG						CELL2022-008				
CHECKED			FINISH: N/A			REFERENCE DRAWING NO.			SCALE: N/A	
APPROVED									MATERIAL:	

PC2022-001



SECTION X - X



Notes:



1. Product data: PTFE, WHITE
2. CTQ (critical to quality) are: Dimension A,B,C,D,E,F and material.
3. Major defects: AQL 1.0 per General Inspection Plan Level II, single sampling (ANSI/ASQC Z1.4-1993)
 - a. Stopper appear dirty or damaged
 - b. Dimension is out of tolerance
4. First Article Report required on new product.
5. Supplier is required to inform VWR and obtain approval before any changes to materials or source, process change or manufacturing location change.
6. Must be manufactured at supplier's premises condition and not from 3rd party.

TABLE OF INFORMATION

Number	DIMENSION(mm)						Cap Mat'l	Suitable for VWR item
	A	B	C	D	E	F		
PC2022-001-1	3.5±0.2	0.8±0.2	5±0.2	2.5±0.2	12.5±0.2	9.6±0.2	PTFE	VWRI634-1196
PC2022-001-2	4.5±0.2	1.8±0.2	5±0.2	2.5±0.2	12.5±0.2	9.6±0.2	PTFE	VWRI634-1197
PC2022-001-5	7.5±0.2	4.6±0.2	5±0.2	2.5±0.2	12.5±0.2	9.6±0.2	PTFE	VWRI634-1198
PC2022-001-10	12.5±0.2	9.6±0.2	5±0.2	2.5±0.2	12.5±0.2	9.6±0.2	PTFE	VWRI634-1199
PC2022-001-10	12.5±0.2	9.6±0.2	5±0.2	2.5±0.2	12.5±0.2	9.6±0.2	PTFE	VWRI634-1200
PC2022-001-20	22.5±0.2	19.6±0.2	5±0.2	2.5±0.2	12.5±0.2	9.6±0.2	PTFE	VWRI634-1201
PC2022-001-40	42.5±0.2	39.6±0.2	5±0.2	2.5±0.2	12.5±0.2	9.6±0.2	PTFE	VWRI634-1202
PC2022-001-50	52.5±0.2	49.6±0.2	5±0.2	2.5±0.2	12.5±0.2	9.6±0.2	PTFE	VWRI634-1203
PC2022-001-100	102.5±0.2	99.6±0.2	5±0.2	2.5±0.2	12.5±0.2	9.6±0.2	PTFE	VWRI634-1204
PC2022-001-10	12.5±0.2	9.6±0.2	5±0.2	2.5±0.2	12.5±0.2	9.6±0.2	PTFE	VWRI634-1205
PC2022-001-10	12.5±0.2	9.6±0.2	5±0.2	2.5±0.2	12.5±0.2	9.6±0.2	PTFE	VWRI634-1206
PC2022-001-20	22.5±0.2	19.6±0.2	5±0.2	2.5±0.2	12.5±0.2	9.6±0.2	PTFE	VWRI634-1207
PC2022-001-40	42.5±0.2	39.6±0.2	5±0.2	2.5±0.2	12.5±0.2	9.6±0.2	PTFE	VWRI634-1208
PC2022-001-50	52.5±0.2	49.6±0.2	5±0.2	2.5±0.2	12.5±0.2	9.6±0.2	PTFE	VWRI634-1209
PC2022-001-100	102.5±0.2	99.6±0.2	5±0.2	2.5±0.2	12.5±0.2	9.6±0.2	PTFE	VWRI634-1210
PC2022-001-10	12.5±0.2	9.6±0.2	5±0.2	2.5±0.2	12.5±0.2	9.6±0.2	PTFE	VWRI634-1216
PC2022-001-20	22.5±0.2	19.6±0.2	5±0.2	2.5±0.2	12.5±0.2	9.6±0.2	PTFE	VWRI634-1217
PC2022-001-50	52.5±0.2	49.6±0.2	5±0.2	2.5±0.2	12.5±0.2	9.6±0.2	PTFE	VWRI634-1218
PC2022-001-10	12.5±0.2	9.6±0.2	5±0.2	2.5±0.2	12.5±0.2	9.6±0.2	PTFE	VWRI634-1219
PC2022-001-10	12.5±0.2	9.6±0.2	5±0.2	2.5±0.2	12.5±0.2	9.6±0.2	PTFE	VWRI634-1223
PC2022-001-20	22.5±0.2	19.6±0.2	5±0.2	2.5±0.2	12.5±0.2	9.6±0.2	PTFE	VWRI634-1224

MARK	ECO NO.		THIS DOCUMENT CONTAINS CONFIDENTIAL AND PROPRIETARY INFORMATION AND MAY NOT BE COPIED OR REPRODUCED IN WHOLE OR IN PART NOR PUT ON CIRCULATIONS WITHOUT THE EXPRESS WRITTEN PERMISSION OF VWR INTERNATIONAL, INC.	ZONE	EDITED	DATE	APPD	DATE	REV
APPROVALS	BY	DATE							
DRAWN		10-Apr-2022	UNLESS OTHERWISE SPECIFIED TOLERANCES AS FOLLOWS DIMENSIONS X ± 0.5 (XX ± 0.25) - DIMENSIONS ARE IN MM UNLESS OTHERWISE SPECIFIED XX ± 0.25 (XXX ± 0.1) - SURFACES FINISH TO BE  XXX ± 0.25 (XXXX ± 0.05) - DIMENSIONS TO BE  ANGLES X ± 1° - WORKMANSHIP SHALL BE THE BEST COMMERCIAL QUALITY XX ± 5°	DESCRIPTION: FLAT PTFE LID					
DESIGN				VWR P/N PC2022-001 REV. A0 SIZE: A3					
MFG				REFERENCE DRAWING NO. SCALE: N/A					
CHECKED				MATERIAL: PTFE, WHITE WEIGHT: 1.0 SHEET 1 OF 1					
APPROVED									

PC2022-002



SECTION X - X




TABLE OF INFORMATION



Number	DIMENSION(mm)					Cap Mat'l	Suitable for VWR item
	A	B	C	D	E		
PC2022-002-1	Ø5±0.2	Ø3.3±0.2	Ø2.5±0.2	8±0.2	5±0.2	PTFE	VWRI634-1211
PC2022-002-1	Ø5±0.2	Ø3.3±0.2	Ø2.5±0.2	8±0.2	5±0.2	PTFE	VWRI634-1212

Notes:

- Product data: PTFE,WHITE
- CTQ (critical to quality) are: Dimension A,B,C,D,E and material.
- Major defects: AQL 1.0 per General Inspection Plan Level II, single sampling (ANSI/ASQC Z1.4-1993)
 - Stopper appear dirty or damaged
 - Dimension is out of tolerance
- First Article Report required on new product.
- Supplier is required to inform VWR and obtain approval before any changes to materials or source, process change or manufacturing location change.
- Must be manufactured at supplier's premises condition and not from 3rd party.

MARK		ECO NO.		THIS DOCUMENT CONTAINS CONFIDENTIAL AND PROPRIETY) INFORMATION AND MAY NOT BE COPIED OR REPRODUCED IN WHOLE OR IN PART NOR PUT ON CIRCULATIONS WITHOUT THE EXPRESS WRITTEN PERMISSION OF VWR INTERNATIONAL, INC	ZONE	EDITED	DATE	APPD	DATE	REV	
APPROVALS	BY	DATE				DESCRIPTION: PTFE STOPPER					
DRAWN		10-Apr-2022		UNLESS OTHERWISE SPECIFIED TOLERANCES AS FOLLOWS DIMENSIONS X ± 0.5 (XXX ± 0.02) XX ± 0.25 (XXX ± 0.01) XXX ± 0.12 (XXX ± 0.005) ANGLES X ± 1° XX ± 0° - DIMENSIONS ARE IN MM UNLESS OTHERWISE SPECIFIED - SURFACES FINISHES TO BE  - DENOTES CHANGE Δ - WORKMANSHIP SHALL BE THE BEST COMMERCIAL QUALITY		VWR P/N			REV. A0		SIZE: A3
DESIGN					REFERENCE DRAWING NO.	PC2022-002					SCALE: N/A
MFG					MATERIAL: PTFE,WHITE				WEIGHT:	SHEET 1 OF 1	
CHECKED					FINISH: N/A						
APPROVED											

PC2022-003



SECTION X - X






TABLE OF INFORMATION

Number	DIMENSION(mm)							Cap Mat'l	Suitable for VWR item
	A	B	C	D	E	F	G		
PC2022-003-5	Ø8±0.2	Ø4±0.2	Ø5±0.2	Ø4.4±0.2	10±0.2	2.7±0.2	6±0.2	PTFE	VWRI634-1213
PC2022-003-10	Ø13±0.2	Ø8±0.2	Ø9.7±0.2	Ø9±0.2	13±0.2	3.5±0.2	7±0.2	PTFE	VWRI634-1214
PC2022-003-10	Ø13±0.2	Ø8±0.2	Ø9.7±0.2	Ø9±0.2	13±0.2	3.5±0.2	7±0.2	PTFE	VWRI634-1215
PC2022-003-10	Ø13±0.2	Ø8±0.2	Ø9.7±0.2	Ø9±0.2	13±0.2	3.5±0.2	7±0.2	PTFE	VWRI634-1225

Notes:

1. Product data: PTFE,WHITE
2. CTQ (critical to quality) are: Dimension A,B,C,D,E,F,G and material.
3. Major defects: AQL 1.0 per General Inspection Plan Level II, single sampling (ANSI/ASQC Z1.4-1993)
4. First Article Report required on new product.
5. Supplier is required to inform VWR and obtain approval before any changes to materials or source, process change or manufacturing location change.
6. Must be manufactured at supplier's premises condition and not from 3rd party.

MARK		ECO NO.		THIS DOCUMENT CONTAINS CONFIDENTIAL AND PROPRIETY) INFORMATION AND MAY NOT BE COPIED OR REPRODUCED IN WHOLE OR IN PART NOR PUT ON CIRCULATIONS WITHOUT THE EXPRESS WRITTEN PERMISSION OF VWR INTERNATIONAL, INC	ZONE	EDITED	DATE	APPD	DATE	REV
APPROVALS	BY	DATE			PROJ: 	DESCRIPTION: PTFE STOPPER				
DRAWN		10-Apr-2022		UNLESS OTHERWISE SPECIFIED TOLERANCES AS FOLLOWS DECIMALS X ± 0.5 (XX ± 0.2) XX ± 0.25 (XXX ± 0.1) XXX ± 0.2 (XXXX ± 0.05) ANGLES X ± 1° XX ± 3° - DIMENSIONS ARE IN MM UNLESS - SURFACES FINISHES TO BE  - DENOTES CHANGE  - WORKMANSHIP SHALL BE THE BEST COMMERCIAL QUALITY	VWR P/N				REV. A0	SIZE: A3
DESIGN					PC2022-003					
MFG					REFERENCE DRAWING NO.	SCALE: N/A				
CHECKED										
APPROVED										
				FINISH: N/A	MATERIAL: PTFE,WHITE	WEIGHT:		SHEET 1 OF 1		



SECTION X - X





TABLE OF INFORMATION


Number	DIMENSION(mm)								Cap Mat'l	Suitable for VWR item
	A	B	C	D	E	F	G	H		
PC2022-004-10	12.5±0.2	9.6±0.2	6.7±0.2	5±0.2	2±0.2	12.5±0.2	9.6±0.2	6.7±0.2	PTFE	VWRI634-1220
PC2022-004-10	12.5±0.2	9.6±0.2	6.7±0.2	5±0.2	2±0.2	12.5±0.2	9.6±0.2	6.7±0.2	PTFE	VWRI634-1221

Notes:

- Product data: PTFE, WHITE
- CTQ (critical to quality) are: Dimension A,B,C,D,E,F,G,H and material.
- Major defects: AQL 1.0 per General Inspection Plan Level II, single sampling (ANSI/ASQC Z1.4-1993)
 - Stopper appear dirty or damaged
 - Dimension is out of tolerance
- First Article Report required on new product.
- Supplier is required to inform VWR and obtain approval before any changes to materials or source, process change or manufacturing location change.
- Must be manufactured at supplier's premises condition and not from 3rd party.

MARK		ECO NO.		THIS DOCUMENT CONTAINS CONFIDENTIAL AND PROPRIETARY INFORMATION AND MAY NOT BE COPIED OR REPRODUCED IN WHOLE OR IN PART NOR PUT ON CIRCULATIONS WITHOUT THE EXPRESS WRITTEN PERMISSION OF VWR INTERNATIONAL, INC.	ZONE	EDITED	DATE	APPD	DATE	REV	
APPROVALS		BY	DATE			DESCRIPTION: FLAT PTFE LID					
DRAWN			10-Apr-2022								
DESIGN				UNLESS OTHERWISE SPECIFIED TOLERANCES AS FOLLOWS DIMENSIONS X ± 0.5 (XXX ± 0.02) XX ± 0.25 (XXX ± 0.01) XXX ± 0.2 (XXX ± 0.05) ANGLES X ± 1° XX ± 5°	VWR P/N PC2022-004 REV. A0 SIZE: A3						
MFG											
CHECKED				- DIMENSIONS ARE IN MM (INCH) - SURFACES FINISHES TO BE  - DENOTES CHANGE  - WORKMANSHIP SHALL BE THE BEST COMMERCIAL QUALITY	REFERENCE DRAWING NO. SHEET N/A						
APPROVED											
				FINISH	MATERIAL		PTFE, WHITE		WEIGHT	SHEET 1 OF 1	

PC2022-005



SECTION X - X




TABLE OF INFORMATION

Number	DIMENSION(mm)							Cap Mat'l	Suitable for VWR item
	A	B	C	D	E	F	G		
PC2022-005-50	52.5±0.2	49.6±0.2	5±0.2	2±0.2	12.5±0.2	9.6±0.2	6.7±0.2	PTFE	VWRI634-1222

Notes:

1. Product data: PTFE,WHITE
2. CTQ (critical to quality) are: Dimension A,B,C,D,E,F,G and material.
3. Major defects: AQL 1.0 per General Inspection Plan Level II, single sampling(ANSI/ASQCZ1.4-1993)
 - a. Stopper appear dirty or damaged
 - b. Dimension is out of tolerance
4. First Article Report required on new product.
5. Supplier is required to inform VWR and obtain approval before any changes to materials or source, process change or manufacturing location change.
6. Must be manufactured at supplier's premises condition and not from 3rd party.

MARK		ECO NO.		THIS DOCUMENT CONTAINS CONFIDENTIAL AND PROPRIETARY INFORMATION AND MAY NOT BE COPIED OR REPRODUCED IN WHOLE OR IN PART NOR PUT ON CIRCULATIONS WITHOUT THE EXPRESS WRITTEN PERMISSION OF VWR INTERNATIONAL, INC.	ZONE	EDITED	DATE	APPD	DATE	REV	
APPROVALS		BY	DATE		PROJ: 						
DRAWN			10-Apr-2022	UNLESS OTHERWISE SPECIFIED TOLERANCES AS FOLLOWS DECIMALS X ± 0.5 XXX ± 0.01 XX ± 0.25 XXX ± 0.12 ANGLES X ± 1° XX ± 5° - DIMENSIONS ARE IN MM UNLESS SURFACES FINISHES TO BE 32 - DIMENTS CHANGE Δ - DIMENSIONS SHALL BE THE BEST COMMERCIAL QUALITY	DESCRIPTION: FLAT PTFE LID						
DESIGN											
MFG					VWR P/N PC2022-005 REV. A0 SIZE: A3						
CHECKED					REFERENCE DRAWING NO. SHEET N/A						
APPROVED				MATERIAL: PTFE,WHITE WEIGHT: 1.0 SHEET 1 OF 1							