

Technical Datasheet

VWR® BASIC PROTECTION DISPOSABLE HOOD

EU cat. No 113-9999

Characteristics


VWR® Basic Protection Disposable non-woven hood is made of 20gsm SPP fabric, and features a small open face hole. Hoods provide reliable protection in laboratory, electrical, pharmaceutical, food, industrial, and cleanroom applications.

- Full garment dimensions are 45.7W x 50.8L cm (18 x 20")
- Constructed from 20gsm SPP fabric

General Specifications

Item	Test method	Standard
Weight	ASTM D3776	20g/m ²
Break Tensile Properties-MD	ASTM D5034	>40 N
Break Tensile Properties-CD	ASTM D5034	>20 N
Elongation Rate-MD	ASTM D5034	>60%
Elongation Rate-CD	ASTM D5034	>70%

Drawing and Dimensions



Item	A	B	C	D	E	F	G	H	I
Dimension (cm)	51	46	21	4	14	35	15	25	15

Packaging and Storage

Package/Quantity	100pcs/bag, 5 bags/carton case
Best storage Conditions	Room temperature, dry and shady
Country of Origin	China

Date of Issue: May 18, 2022