

Certificate of Conformity

VWR Catalogue Number	See List
Description	Carboys, narrow neck, round, with stopcock and screw cap
Country of Origin	Manufactured in India
Date of Issue (yyyy-mm-dd)	2022-05-23

Low Density Polyethylene Carboys

Catalogue Number	Capacity [L]	Thread [mm]
216-1909	5	53
216-1731	10	83
216-1717	20	83
216-1899	25	83
216-1734	50	83

High Density Polyethylene Carboys

Catalogue Number	Capacity [L]	Thread [mm]
216-1928	25	83

Polypropylene Carboys

Catalogue Number	Capacity [L]	Thread [mm]
216-1710	10	83
216-1705	20	83
216-1718	50	83

Replacement caps

Catalogue Number	Type	Cap Size [mm]	Colour
216-1840	Cap for narrow mouth VWR® carboys	83	Purple

Replacement stopcocks

Catalogue Number	Type
216-1745	Stopcock, with adapter for tubing connection and stopper for long storage for VWR Collection carboys

Quality System Compliance

Products are manufactured under the **ISO 9001:2015 & ISO 13485:2016** standard. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements.

For LDPE and HDPE Carboys

Test Name	Test Reference	Accuracy	Result
Leak Test	QWI-QC-04	General Insp. Level-1, AQL -1.5	Leak Proof
Autoclaving	QWI-QC-07	General Insp. Level-1, AQL -1.5	NA
Dimensional	QSF-QCIR-03	General Insp. Level-1, AQL -1.5	Within Tolerance

For PP Carboy

Test Name	Test Reference	Accuracy	Result
Leak Test	QWI-QC-04	General Insp. Level-1, AQL -1.5	Leak Proof
Autoclaving	QWI-QC-07	General Insp. Level-1, AQL -1.5	PASS
Dimensional	QSF-QCIR-03	General Insp. Level-1, AQL -1.5	Within Tolerance

For Replacement caps

Test Name	Test Reference	Accuracy	Result
Leak Test	QWI-QC-04	General Insp. Level-1, AQL -1.5	Leak Proof
Autoclaving	QWI-QC-07	General Insp. Level-1, AQL -1.5	PASS
Dimensional	QWI-QC-23	General Insp. Level-1, AQL -1.5	Within Tolerance

For Replacement stopcocks

Test Name	Test Reference	Accuracy	Result
Leak Test	QWI-QC-19	General Insp. Level-1, AQL -1.5	Leak Proof
Autoclaving	QWI-QC-07	General Insp. Level-1, AQL -1.5	PASS
Dimensional	QWI-QC-23	General Insp. Level-1, AQL -1.5	Within Tolerance

Product Specifications

Material	Carboy: Low Density Polyethylene or Polypropylene, 100 % Virgin, FDA 21CFR Compliance, USP Class VI passed. Closure material: Polypropylene, 100 % Virgin, FDA 21CFR Compliance, USP Class VI passed with Food grade colour.	
Temperature Resistance	Low Density Polyethylene -80 to +80 °C. High Density Polyethylene -80 to +120 °C. Polypropylene 0 to +135 °C.	
Sterilization	Not Applicable.	
ATP Assay	Not Applicable.	
DNase & RNase Free	Not Applicable.	
BSE/TSE	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).	
Cytotoxicity	Not Applicable.	
Non-Pyrogenic Statement	Not Applicable.	
Latex Statement	This product is not made from natural rubber latex.	
BPA Statement	Not Applicable.	
DEHP Statement	Please see RoHS information below.	

RoHS	No substances (Lead, Cadmium, Mercury, Hexavalent Chromium (Cr6+), Poly Brominated Biphenyls (PBB), Poly Brominated Diphenyl ethers (PBDE), Bis(2-Ethyhexyl) Phthalate (DEHP), Benzyl Butyl Phthalate (BBP), Dibutyl Phthalate (DBP), Diisobutyl Phthalate (DIBP)) are used in manufacturing the raw materials and final product. No routinely analyse is performed.
REACH Statement	No substances as listed on >Substance of Very High Concern< candidate list published 2022-01-17 are used in manufacturing of the raw materials used to produce this product. No routinely analyse is performed.
Storage Conditions	Ambient temperature.
Shelf Life	Not Applicable.

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
This document has been produced electronically and is valid without a signature.