



Certificate of Conformance

Sharps Container

Item Number:	19001-001
Lot Numbers:	494543
MFG Date	7/3/2025
Description	VWR Container System 1 QT
Material:	Polypropylene
Sterile:	No
Expiration Date:	N/A
Country of Origin:	USA

Non Pyrogenic

Product has not been treated for pyrogens

DNase & RNase Free

Product has not been treated for DNase or RNase agents

Cytotoxicity

N/A

BSE/TSE: Resins used to mold the sharps container may contain additives from North American cattle. Any tallow used has been treated by transesterification or hydrolysis at 200 degrees C for 20 minutes at 40 bar or greater.

This is to certify that this product has been inspected and released in accordance with Medegen Medical Products Standard Production Methods and Process Quality Assurance Specifications. All resins used in the processing of these products have been approved and released. The product was manufactured in compliance with current cGMP's.

Authorized Signature: Raven Wade

Date: 9/2/2025