

Certificate of Conformity

VWR Catalogue Number	See List
Description	Lab coats, Basic, SMS, disposable, VWR®
Country of Origin	Manufactured in China
Date of Issue (yyyy-mm-dd)	2025-04-11

Catalogue Number	Colour	Clothing size	Length [cm]
113-1131	White	S	93
113-1132	White	M	100
113-1133	White	L	102
113-1134	White	XL	108
113-1135	White	XXL	110
113-1136	White	3XL	113
113-1137	White	4XL	118
113-1125	Blue	S	93
113-1126	Blue	M	100
113-1127	Blue	L	102
113-1128	Blue	XL	108
113-1129	Blue	XXL	110
113-1130	Blue	3XL	113

Quality System Compliance

Products are manufactured under the **ISO 9001:2015 & ISO 13485:2016** standard. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements. Inspected according to AQL level II 4.0 sampling plan per MIL-STD-105E.

	Test Method	MD	CD
Break Tensile Properties [N/5cm]	ASTM D5034 standard	> 70	> 60
Elongation rate [%]	ASTM D5034 standard	> 50	> 50

Product Specifications

Material	SPP 50 g/m ² (Spunbonded Polypropylene fabric+Melt bown+Spunbonded Polypropylene fabric).
Sterilization	Non-sterile.
ATP Assay	Not Applicable.
DNase & RNase Free	Not Applicable.
BSE/TSE	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
Cytotoxicity	Not Applicable.
Non-Pyrogenic Statement	Not Applicable.
Latex Statement	This product is not made from natural rubber latex.
BPA Statement	Not Applicable.
DEHP Statement	Not Applicable.
RoHS	Not Applicable.
REACH Statement	Not Applicable.
Storage Conditions	Room temperature, dry and shady.
Shelf Life	5 years.

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
This document has been produced electronically and is valid without a signature.

Version History

Sl.no	Version number	Description	Month & Year
1	1.0	Created General CoQ using new layout	11-2021
2	1.1	Updated Shelf-life information	04-2025