

Calbiochem[®]

Certificate of Analysis


Leptin, Human, Recombinant, *E. coli*

Batch Number:	D00024276
Material Number:	429700-1MG
Molecular Weight:	16000
CAS Number:	177404-21-6

Analytical Data

Test	Tolerance	Result
Formulation:		Lyophilized from a sterile-filtered solution in PBS
Form:		Solid
Purity by SDS-PAGE:	≥97.0 %	>97.0 %
Endotoxins:	≤1.0 EU/μg	<1.0 EU/μg

Storage and Handling:	≤ -70C
------------------------------	--------



Ellen Kroninger, Quality Assurance Manager

19-May-2008

Date

USA and Canada
Tel (800) 628-8470
technical@calbiochem.com

France
Freephone
0800 126 461

Germany
Freecall
0800 100 3496

Ireland
Toll Free
1800 409 445

Europe
United Kingdom
Freephone
0800 622 935

All Other
European Countries
+44 115 943 0840

All Other Countries
Contact Your Local
Distributor
www.calbiochem.com
technical@calbiochem.com

techservice@merckbio.eu

A Brand of EMD Chemicals Inc., an affiliate of Merck KGaA, Darmstadt, Germany
www.calbiochem.com

FOR RESEARCH USE ONLY. NOT FOR HUMAN OR DIAGNOSTIC USE.