

Data sheet

Product: ringcaps® 5 + 10 µl

Code-no.: 960 02 10

The above named product conforms with the requirements of the standard ISO 7550.

Glass material: DURAN, borosilicate glass 3.3 (acc. ISO 3585)

Inner diameter: $0,47 \pm 0,05$ mm

Outer diameter: $1,20 \pm 0,05$ mm

Tolerances: Accuracy A: $\pm 0,25\%$
Coefficient of variation: $0,50\%$

Graduation line: 5 + 10 µl

Color code: Orange

Eberstadt, 30.01.2023



.....
Roland Benz/ Quality Manager