

Technical Data Sheet



VWR® Centrifuge Tube

Description

50mL, PS, Conical Bottom, Flat Cap, Sterilized

Purpose

Useful for collecting and transporting biological specimens

Materials

Tube: PS (Polystyrene)

Color: Clear




Cap: HDPE (High-density Polyethylene)

Color: White

Rack: PP (Polypropylene)

Color: White

Dimensions (Unit: mm)



Features

- Made with high-quality USP Class VI polystyrene
- Polystyrene tubes can tolerate aqueous solutions of mild bases or weak acids, but not organic solvents, aromatic or chlorinated hydrocarbons
- Polystyrene tubes are not autoclavable
- Maximum RCF up to 2,000 x g
- Rack marked with numbers (1~25) for identification
- Rack can be washed and reused
- DNase/RNase free
- Non-pyrogenic
- Sterilized by irradiation SAL 10⁻⁶ (ISO11137)
- Shelf Life: 3 years after month of production
- Manufactured in a class 100,000 clean room environment
- Manufactured under ISO13485 and ISO9001 quality management system



VWR® Centrifuge Tube

VWR NA Cat. No.	Capacity	Maximum RCF	Bottom Style	Cap Style	Sterilization	Packaging Configuration
78592-328	50 mL	2,000 x g	Conical Bottom	Flat cap	Yes	25/Rack/Bag, 300/CS