

Certificate of Conformity

Catalog Number	07-5050
Product Description	Serological Pipette 50ml, purple plug, Sterile to SAL10 ⁻⁶ , 1/bag, 40bags/pack, 4bags/case
Lot Number	CB1E09A
Expiration Date	May., 2026

Visual Attribute	Pass
------------------	------

Functional Test	Pass
-----------------	------

Sterilization	Sterilization was performed by exposed electron beam according to ISO 11137 based on the Sterilization of healthcare products-Requirements for validation and routine control-Radiation sterilization <SAL 10 ⁻⁶ >.
---------------	--

Cytotoxicity Test	Cytotoxicity test is conducted to qualify all material resins using USP class VI <87> and/or ISO 10993 standards for cytotoxicity and have been shown to be non-toxic. endotoxin-free
-------------------	---

Pyrogens Test	Non-pyrogenic tested according to the principles of the LAL test described in USP <85> / 0.01EU/ml.
---------------	---

RNase/DNase/DNA Test	This product has been tested and is free of any detectable RNase/DNase/DNA contamination.
----------------------	---

Biologix Group Ltd. hereby certifies that the product identified above have been inspected in compliance with product quality specification and requirements as documented in the Quality System Standard ISO 9001:2015.

Date:	02 November 2023
-------	------------------

Authorized Signature

