

Certificate of Conformity

VWR Catalogue Number	612-1701
Description	VWR®, Pasteur Pipettes, non plugged, 2 mL capacity, 150 mm length
Lot Number	180745
Country of Origin	Manufactured in Germany
Date of Issue (yyyy-mm-dd)	2026-04-21

Quality System Compliance

This certificate confirms that the product has been sourced from a supplier or manufacturer who has declared their Quality Management System (QMS) complies with [ISO 9001:2015](#).

QC Testing

Representative products samples are collected and inspected in accordance with current applicable product specifications.

Product Specifications

Top body outer diameter: 7.1 mm

Wall thickness: 0.53 mm

Jet outer diameter: 1.5 mm

Material Soda Lime glass tubing.

Sterilization Non-Sterile.

ATP Assay Not Applicable.

DNase & RNase Free Not Applicable.

BSE/TSE Not Applicable.

Cytotoxicity Not Applicable.

Non-Pyrogenic Statement Not Applicable.

Latex Statement	Not Applicable.
BPA Statement	Not Applicable.
DEHP Statement	Not Applicable.
RoHS	Not Applicable.
REACH Statement	Not Applicable.
Storage Conditions	Storage at ambient temperature.
Shelf Life	Not Applicable.

Disclaimer: VWR's Quality Management System (QMS) is ISO 9001:2015 certified as a wholesale distributor of finished products, ensuring effective oversight of finished product handling, traceability, and distribution processes. VWR is not a manufacturer and does not engage in manufacturing, co-manufacturing or processing activities such as design, packaging, labeling, testing or quality release. Consequently, the supplier or manufacturer's QMS certification may differ from VWR's.

VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. This product is intended for research use only. It is not for use in diagnostic procedures, therapeutic applications, or any human or clinical use.

This document has been produced electronically and is valid without a signature.