



Certificate of Conformity

VWR Conductivity Meter, Handheld, HCO 304

VWR International certifies that the products comply in every respect with the relevant particular specifications and/or drawings and/or relevant standards and regulations in force, mentioned in this document.

VWR International certifie que les produits sont conformes à tous égards aux spécifications particulières appropriées et / ou dessins, normes et règlements en vigueur, mentionnées dans le présent document

VWR International bestätigt, dass das Produkte allen zutreffenden Spezifikations- und/oder Zeichnungs- und/oder Standardanforderungen sowie den im Dokument aufgeführten gültigen Vorschriften entspricht.

Cat N°: 663-0348

Specifications:

Conductivity range: 0 - 2000 $\mu\text{S}/\text{cm}$
0.00 ... 20.00 mS/cm
0.0 ... 100.0 mS/cm

Conductivity accuracy: $\pm 0.5\%$ of m.v. $\pm 0.5\%$ FS

Salinity range: 0.0 – 50.0 g/kg water

TDS range: 0 - 2000 mg/l

Temperature range (°C): -5.0...+100.0

Temperature accuracy (°C): ± 0.3

Temperature compensation : Off: deactivated
nLF: non linear, according to EN 27888

Display: 3-line display with backlight.

Power supply: 2× AA battery

W×D×H (mm): 54×28×108

Weight (kg): 0.2

Compliance: IP 65/IP 67

Sensor type: 2-pole graphite measuring cell

Certificate of Conformity

Other Characteristics:

- Measures conductivity, salinity and TDS.
- Permanently connected 2-pole measuring cell (graphite) with integrated temperature sensor.
- Dust- and waterproof (IP65/IP67) ABS housing.
- 3-line display with battery status indicator and backlight, protected by an unbreakable pane, overhead display at the push of a button.
- HOLD function, min/max values, automatic measuring range shifting, automatic temperature compensation.

Disclaimer: VWR wants to state that this declaration will not discharge the user from his obligation to ensure the product is suitable for the intended use. - If there is any doubt ask us.

Date of issue: **July 2019**