



ACROS ORGANICS part of Thermo Fisher Scientific

ACRÖS
ORGANICS

Product Specifications

Product 131780250
Pyridine, 99+%, extra pure

General Product Data

Version	6
CAS No	110-86-1
Molecular weight	79.10
Molecular formula	C ₅ H ₅ N
Linear formula	
Flash point (°C)	17

Product Specifications

Appearance (Form)	Clear liquid
Color scale	=<40 APHA
Infrared spectrum	Conforms
GC	>=99.0 %
Water	=<0.1 % (Coulometric)
Residue after evaporation	=<0.02 %
Refractive index	1.5080 to 1.5110 (20°C, 589 nm)

Product Specifications are subject to amendment and may change over time.

Acros Organics

ENA23, zone 1, nr 1350, Janssen Pharmaceutica laan 3a, B-2440 Geel, Belgium
Tel +32 14/57.52.11 - Fax +32 14/59.34.34 Internet: <http://www.acros.com>
1 Reagent Lane, Fair Lawn, NJ 07410, USA
Fax 201-796-1329