

Certificate of Quality

VWR Catalogue Number	See List
Description	VWR® Stainless Steel Dissecting Forceps with Fine, Serrated Tips
Country of Origin	Made in Pakistan
Date of Issue (yyyy-mm-dd)	2020-08-20

Catalogue Number	Length	Description
82027-408	114 mm (4 1/2")	Straight
82027-410	114 mm (4 1/2")	Curved

Quality System Compliance

Products are manufactured under the [ISO 9001:2008 & ISO13485:2003 & F-899](#) standard. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements.

Product Specifications

Material	AISI 410 Type Stainless steel fabrication
Sterilization	Not applicable
ATP Assay	Not applicable
DNase & RNase Free	Not applicable
BSE/TSE	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).

Cytotoxicity	Not applicable
Non-Pyrogenic Statement	Not applicable
Latex Statement	The product is Latex free
BPA Statement	Not applicable
DEHP Statement	Not applicable
RoHS	Not applicable
REACH Statement	Not applicable
Storage Conditions	Store at dry place
Shelf Life	Not applicable

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
This document has been produced electronically and is valid without a signature.