



Product Specification Sheet

VWR TEMPERATURE STRIPS, 8-LEVEL

VWR International certifies that the specified products meet all VWR specifications and quality criteria for release. The products comply in every respect with the relevant particular specifications, drawings, relevant standards and regulations in force.

VWR International certifie que les produits mentionnés respectent toutes les spécifications et critères de qualité VWR. Les produits satisfont aux spécifications particulières, aux plans, aux normes et règlements en vigueur s'y rapportant.

VWR International bestätigt hiermit, daß die angegebenen Produkte alle VWR Spezifikationen und Qualitätskriterien zur Freigabe erfüllen. Die Produkte entsprechen in jeder Hinsicht den maßgeblichen Spezifikationen, den Zeichnungen sowie den relevanten Normen und gültigen Vorschriften.

Cat. N°	Temperature Range (°C)	Pack (pces)
620-9102	+37 ... +65	10
620-9103	+71 ... +110	10
620-9104	+116 ... +154	10
620-9105	+160 ... +199	10
620-9106	+204 ... +260	10

Product Description: A series of temperature test labels that will permanently, record, the highest temperature reached by the label.

Application Data:

- Peel indicator from release paper
- Apply to any clean, grease free surface ensuring that the indicator has adhered to the surface
- A color change from silver white/orange to black indicates that the temperature has been reached



Product Specification Sheet

VWR TEMPERATURE STRIPS, 8-LEVEL

Physical Properties:

	Low Temperatures (171°C and below)	Mid Temperatures (177°C and below)	High Temperatures 210°C +
Adhesive Type	Acrylic	Modified Acrylic	Modified Acrylic
Carrier	Polyester	Tissue	Tissue
Covering Film	Polyester	Kaladex	Polyimide
Colour Change Material	Non toxic, white crystalline solid on a black absorbent backing, adhered with acrylic adhesive		
Tolerance	Less than 100°C = $\pm 1^{\circ}\text{C}$, 100°C to 154°C = $\pm 1.5^{\circ}\text{C}$ Greater than 154°C = $\pm 4^{\circ}\text{C}$		
Shelf Life	12 months from invoice date when stored in a cool dry environment (<64°F & < 50% relative humidity)		
Thickness	At temperature indicator chemical 6.5 mil No temperature indicator chemical 1.5mil		

Disclaimer: VWR wants to state that these specifications can only be regarded as general information and will not discharge the user from his obligation to ensure the product is suitable for the intended use - If there's any doubt ask us.