



Product Information

Peroxide-Free Detergents

For proteomics applications.

Tween 20	M228-10ML-5PK
Triton X-100	M236-10ML-5PK
Triton X-114	M237-10ML-5PK
Brij 35	M238-10ML-5PK
Detergent Sampler	M248-4PK

AMRESCO offers a variety of peroxide-free nonionic detergents in a convenient ready-to-use 10% solution. These detergents are purified to remove peroxide contamination (less than 1 ueq/ml peroxides) and then packaged in sealed glass ampoules under argon. Each of the detergents are sold as a 5-pack of 10 ml ampoules. Not sure which detergent works best for your application? AMRESCO also offers a sampler 4-pack with one ampoule of each detergent.

In addition to peroxide-free detergents, AMRESCO offers an extensive line of proteomics reagents and kits. Some of our more popular items include:

Silver BULLit silver stain kit	M229-1L-KIT
BCA protein determination kit	M227-1L-KIT
Proteo-Mix gel kit	M162-KIT, M163-KIT, M164-KIT
Proteo-Mix PLUS gel kit	M175-KIT, M176-KIT, M177-KIT
Protease inhibitor cocktails	M221, M222

ZY0274 06/02

©Copyright 2002 by AMRESCO Inc. All Rights Reserved. AMRESCO® is a registered trademark of AMRESCO Inc.

