



Certificate of Conformity

VWR HAND TALLY COUNTER

VWR International certifies that the products comply in every respect with the relevant particular specifications and/or drawings and/or relevant standards and regulations in force, mentioned in this document.

VWR International certifie que les produits sont conformes à tous égards aux spécifications particulières appropriées et / ou dessins, normes et règlements en vigueur, mentionnées dans le présent document

VWR International bestätigt, dass das Produkte allen zutreffenden Spezifikations- und/oder Zeichnungs- und/oder Standardanforderungen sowie den im Dokument aufgeführten gültigen Vorschriften entspricht.

Cat No: 710-0935

Size: 1-3/4-inch diameter by 1-1/8-inch width

Pk: 1

Other Characteristics:

Compact unit is ideal for carrying in the pocket. Easy-to-use

Four-digit unit records up to 9999. Numbers instantly and automatically return to zero with a twist of the reset knob.

Metal construction is designed for reliability, long life, and hard use. Unique engineering provides a smooth, tactile feel.

Hundreds of laboratory applications provide opportunities to, for example, count blood samples, take surveys, inventory items, total drops of liquid, receive lab samples, and tally biomedical events.

For ease in viewing, large white numerals appear on a black background. A protective window seals out dust, dirt, and humidity.

The rugged, heavy-duty chrome-finished case has a hinged finger loop for a sure grip. Finger pressure on the button produces an audible click that signals the count. This feature provides the assurance of a registered count without looking at the counter. Corrosion-proof stainless-steel for long-life.

Disclaimer: VWR wants to state that this declaration will not discharge the user from his obligation to ensure the product is suitable for the intended use. - If there is any doubt ask us.

Date of issue: **May 2017**