



L3 2063060

OKC Allergy Supplies, Inc.

a division of

ALK-Abello, Inc.

Certificate of Analysis

Product: SEV - Sterile Empty Vial
Lot number: SEV2100521
Vial size: 100mL-20mm
Manufacture date: 10/5/2021
Part Number: SEV100

Quality Release Testing		
Assay	Specifications	Results
Appearance - Product	Essentially free of visible particulate	Pass
Appearance - Packaging	No apparent leakage or physical alteration	Pass
Sterility, USP <71>	Sterile	Conforms
Bacterial Endotoxins, USP <85>	≤ 0.25 EU/mL	≤ 0.0050 EU/mL

Release Date	Expiration Date
10/28/2021	10 - 2025

Quality Assurance Approval	
By / Date: 	11-01-2021

Original Certificate