



The BASICBOX storage boxes are specially designed to meet the needs of laboratories for storing and organizing various types of products. Made of polypropylene, they offer excellent shock resistance while remaining lightweight, making them ideal for intensive use in laboratory settings.

Thanks to their transparency, BASICBOX containers allow for easy identification of the contents, facilitating organization and quick access to samples and reagents. Additionally, they can be cleaned by hand or in a dishwasher, ensuring optimal hygiene, which is essential in a laboratory environment.

The snap-on lid closes securely, protecting the contents from external contaminants. These boxes are designed to be stackable with excellent stacking stability, allowing them to be safely and stably superimposed, thus optimizing storage space.

## PRODUCT ADVANTAGES

- **Secure** : The lid of the box provides constant protection for the contents, safeguarding products from dust, moisture, and other external contaminants.
- **Practical** : Transparent, it is suitable for storing various types of products, making it versatile for multiple applications, whether for storage or transport.
- **Durable** : Made of polypropylene, it is resistant to temperature fluctuations, moisture, and impacts.
- **Hygienic** : Full, smooth walls allow for quick drainage of liquids.

## SPECIFICATIONS

|                       |                   |
|-----------------------|-------------------|
| SKU                   | MULC2291          |
| Ext. dim.             | 205 x 170 x 75 mm |
| Int. dim. (Top)       | 180 x 150 mm      |
| Int. dim. (Bottom)    | 170 x 140 mm      |
| Weight                | 100 g             |
| Volume                | 1.7 L             |
| Material              | Polypropylene     |
| Type of walls / base  | Solid / Solid     |
| Color                 | Transparent       |
| Food Safe             | Yes               |
| Stackable / nestable  | Yes               |
| Recyclable            | Yes               |
| Operating temperature | -20°C / +100°C    |
| Shipping time         | 5 working days    |
| Origin                | Denmark           |