



## Certificate of Conformity

---

### **VWR DIGITAL TIMER, ONE CHANNEL, WHITE**

---

VWR International certifies that the products specified conform to all specifications and quality criteria mentioned in this document.

*VWR International certifie que les produits mentionnés respectent toutes les spécifications et critères de qualité mentionnés dans ce document.*

VWR International bestätigt hiermit, daß die angegebenen Produkte alle in diesem Dokument erwähnten Spezifikationen und Qualitätskriterien erfüllen.

VWR Cat. No: **609-0184**

WxDxH: **46×12×88 mm**

Weight: **35.4 g**

Other Specifications:

- **Easy-to-use timer has a repeat memory function for repetitive tests.**
- **Three-key operation can be learned in less than 10 seconds.**
- **May be placed on incubators, refrigerators, and ovens with the magnet on the back.**
- **Counts down from 99 minutes, 59 seconds in 1-second increments.**
- **Memory retains last programmed time, eliminating the need to re-enter test times.**
- **Alarm may be silenced manually or turns off automatically after 30 seconds.**
- **Supplied with silver-oxide battery, magnet, and Velcro®.**

**Disclaimer:** VWR wants to state that this declaration will not discharge the user from his obligation to ensure the product is suitable for the intended use. - If there is any doubt ask us.

Date of issue: **April 2015**