

Product Specification

Test	Specification
Assay	$\geq 85 \%$
Aluminium (Al)	$\leq 0.001 \%$
Calcium (Ca)	$\leq 0.001 \%$
Potassium Carbonate (K_2CO_3)	$\leq 1.0 \%$
Silica (SiO_2)	$\leq 0.005 \%$
Sodium (Na)	$\leq 0.5 \%$
Trace Impurities (in ppm)	
Chloride (Cl)	≤ 5
Heavy Metals (as Pb)	≤ 5
Iron (Fe)	≤ 5
Nickel (Ni)	≤ 5
Nitrogen Compounds (as N)	≤ 3
Phosphate (PO_4)	≤ 5
Sulphate (SO_4)	≤ 5

For Laboratory, Research or Manufacturing use.



Barbara Rykowska
Quality Assurance Coordinator

For questions on this Product Specification please contact Technical Service at technical.service@avantorsciences.com
Avantor Performance Materials Poland S.A. Sowinskiego 11 Gliwice 44-101 Poland Tel: +48 322392000