

Certificate of Conformity

VWR Burettes, Straight Stopcock

VWR International certifies that the specified products meet all VWR specifications and quality criteria for release. The products comply in every respect with the relevant particular specifications, drawings, relevant standards and regulations in force.

VWR International certifie que les produits mentionnés respectent toutes les spécifications et critères de qualité VWR. Les produits satisfont aux spécifications particulières, aux plans, aux normes et règlements en vigueur s'y rapportant.

VWR International bestätigt hiermit, daß die angegebenen Produkte alle VWR Spezifikationen und Qualitätskriterien zur Freigabe erfüllen. Die Produkte entsprechen in jeder Hinsicht den maßgeblichen Spezifikationen, den Zeichnungen sowie den relevanten Normen und gültigen Vorschriften.

Cat. No.	Lot No.	Capacity ml	Division ml	Tolerance	Scale Length mm	Pack
612-4139	GL2325	10	0.05	0.02 $\pm$ ml	500	1

Date of Manufacture: 02nd June 2025

Quality System Compliance: This certificate confirms that the product has been sourced from a supplier or manufacturer who has declared their Quality Management System (QMS) complies with ISO 9001:2015 standard and calibrated in ISO 17025:2017.

QC Testing: Representative production samples are collected and inspected in accordance with current applicable product specifications.

Product Specifications:

The product is in conformance with the requirements of the following standards and regulations:

Technical standards for products: EN ISO 385

Calibration procedure: ISO 4787

Material:

The products are made from DIN ISO 3585:1998, TYPE-1, CLASS-A, Boro3.3 Glass that meet ASTM E-438 TYPE-1, CLASS-A standard.

BPA statement: Not Applicable

Phthalate/DEHP Statement: Not Applicable

Heavy metals-RoHS: Not Applicable

Certificate of Conformity

REACH (SVHC) Statement: Not Applicable

BSE/TSE statement: The product is solely manufactured from Glass and is free from BSE and TSE. Glass is manufactured by using raw material of non-animal origin and does not contain meat product and does not come in contact with meat products.

Non-Pyrogenic: Not Applicable

DNase & RNase Free: Not Applicable

Cytotoxicity: Not Applicable

Latex Free statement: Not Applicable

Storage Conditions: Can be stored at room temperature in normal warehouse conditions

Shelf Life: There is no predefined shelf life assigned to the product.

Additional Information: The product is safe to use when used in proper way, as should be use with lab glassware. The purpose of the product is for use in laboratory only.

Disclaimer: VWR's Quality Management System (QMS) is ISO 9001:2015 certified as a wholesale distributor of finished products, ensuring effective oversight of finished product handling, traceability, and distribution processes. VWR is not a manufacture and does not engage in manufacturing, co-manufacturing or processing activities such as design, packaging, labeling, testing or quality release. Consequently, the supplier or manufacturer's QMS certification may differ from VWR's.

VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. This product is intended for research use only. It is not for use in diagnostic procedures, therapeutic applications, or any human or clinical use.

This document has been produced electronically and is valid without a signature.

Date of issue: **July 2025**