

CERTIFICATE OF CONFORMITY

PART #:	LOT #:
ASWF1S	20240504A-12A
DESCRIPTION:	
SMALL DISPOSABLE ANTI-STATIC POLYPROPYLENE WEIGING FUNNEL. FITS INTO 5-50mL Volumetric Flasks. *INTENDED FOR SINGLE USE ONLY*	
SPECIFICATIONS:	MATERIAL:
LENGTH: 43.18mm FUNNEL O.D.: 22.86mm POURING NECK O.D.: 7.62mm CAPACITY: 2mL SOLID/2mL LIQUID	VIRGIN POLYPROPYLENE

TWD Scientific, LLC certifies that:

1. The product will not be altered or modified in any manner from the original condition after the product(s) are manufactured and packaged.
2. The product(s) meet the description specifications as detailed in the product description above.

Prepared by: Christina Weber / Customer Service

