



AquaPhoenix®
S C I E N T I F I C



CHEMetrics®
By AquaPhoenix Scientific

Simplicity in Water Analysis

Certificate of Conformance

Catalog No: K-5510

Product Description: Peroxide / Persulfate Kit

Lot: 4A051

Date of Manufacture: 27 February 2024

Date Certified: 27 February 2024

Expiration Date: N/A

Test	Specification	Answer
Reagent Appearance	colorless	Pass
Hydrogen Peroxide Test at 0 ppm	0 to 0.1 ppm	Pass
Hydrogen Peroxide Test at 8.0 ppm	7 to 10 ppm	Pass

CHEMetrics by AquaPhoenix Scientific, LLC hereby certifies that the following water analysis products have been manufactured and evaluated in accordance with strict Quality Control (QC) and Quality Assurance (QA) protocols and/or conform to manufacturer specifications. Where applicable, laboratory test results as well as retention samples for each production batch are on file at CHEMetrics.

Jennifer R. Nash

Signature

Jennifer R. Nash
Director of Quality



AquaPhoenix[®]
S C I E N T I F I C



CHEMetrics[®]
By AquaPhoenix Scientific

Simplicity in Water Analysis

Certificate of Conformance

Catalog No: WC-5501

Product Description: Hydrogen Peroxide Comp. (round)

Lot: 3GI1360

Date of Manufacture: 02 October 2023

Date Certified: 02 October 2023

Expiration Date: September 2026

Test	Specification	Answer
Reagent Appearance	gradient from colorless to light orange then to dai	PASS

CHEMetrics by AquaPhoenix Scientific, LLC hereby certifies that the following water analysis products have been manufactured and evaluated in accordance with strict Quality Control (QC) and Quality Assurance (QA) protocols and/or conform to manufacturer specifications. Where applicable, laboratory test results as well as retention samples for each production batch are on file at CHEMetrics.

Jennifer R. Nash

Signature

Jennifer R. Nash
Director of Quality



AquaPhoenix[®]
S C I E N T I F I C



CHEMetrics[®]
By AquaPhoenix Scientific

Simplicity in Water Analysis

Certificate of Conformance

Catalog No: WC-5510

Product Description: Hydrogen Peroxide Comp. (flat)

Lot: 3GI1360

Date of Manufacture: 02 October 2023

Date Certified: 02 October 2023

Expiration Date: September 2026

Test	Specification	Answer
Reagent Appearance	gradient from light orange to dark orange	PASS

CHEMetrics by AquaPhoenix Scientific, LLC hereby certifies that the following water analysis products have been manufactured and evaluated in accordance with strict Quality Control (QC) and Quality Assurance (QA) protocols and/or conform to manufacturer specifications. Where applicable, laboratory test results as well as retention samples for each production batch are on file at CHEMetrics.

Jennifer R. Nash

Signature

Jennifer R. Nash
Director of Quality



AquaPhoenix®
S C I E N T I F I C



CHEMetrics®
By AquaPhoenix Scientific

Simplicity in Water Analysis

Certificate of Conformance

Catalog No: WR-5510

Product Description: Peroxide / Persulfate

Lot: 4GA1231

Date of Manufacture: 01 February 2024

Date Certified: 01 February 2024

Expiration Date: July 2026

Test	Specification	Answer
Reagent Appearance	colorless	Pass
Hydrogen Peroxide Test at 0 ppm	0 to 0.1 ppm	Pass
Hydrogen Peroxide Test at 8.0 ppm	7 to 10 ppm	Pass

CHEMetrics by AquaPhoenix Scientific, LLC hereby certifies that the following water analysis products have been manufactured and evaluated in accordance with strict Quality Control (QC) and Quality Assurance (QA) protocols and/or conform to manufacturer specifications. Where applicable, laboratory test results as well as retention samples for each production batch are on file at CHEMetrics.

Jennifer R. Nash

Signature

Jennifer R. Nash
Director of Quality