

Product Specification

Test	Specification
Styrene Divinylbenzene (SDB-1):	
Breakthrough Antrazine*	≥ 1 L
Breakthrough Desisopropylatrazine*	≥ 1 L
Specific Surface Area	≥ 950 m ² /g
Porosity, Median Pore Diameter	50 – 70 Å
Particle Size – 50% Volume	65 – 95 µm
Capacity (mg caffeine / g sorbent)	≥ 250
Recovery Atrazine*	≥ 97 %
Recovery Desisopropylatrazine*	≥ 75 %

For Laboratory, Research or Manufacturing use.
1mL Solid Phase Extraction Columns
50 mg per Column



Barbara Rykowska
Quality Assurance Coordinator

For questions on this Product Specification please contact Technical Service at technical.service@avantorsciences.com
Avantor Performance Materials Poland S.A. Sowinskiego 11 Gliwice 44-101 Poland Tel: +48 322392000