

Certificate of Quality

VWR Catalogue Number	76371-722
Description	VWR® TAE Buffer Powder Packs
Country of Origin	Manufactured in China
Date of Issue (yyyy-mm-dd)	2021-08-19

Quality System Compliance

Products are manufactured in accordance with the manufacture's quality system. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable product specifications.

Assay Name	Criteria	Result
Weight (g/pouch)	7.12 (7.82 - 7.42)	Pass
Appearance	White granules	Pass
pH at 25°C	8.2 - 8.4	Pass
Dissolving time at 25°C	Less than 5 minutes	Pass

Product Specifications

Material	40mM Tris, 20mM acetic acid, 1mM EDTA.
Sterilization	Not applicable.
ATP Assay	Not applicable.
DNase & RNase Free	This product is free of any detectable DNase/RNase contamination.
Protease	This product is free of any protease contamination
BSE/TSE	Not applicable.
Cytotoxicity	Not applicable.
Non-Pyrogenic Statement	Not applicable.
Latex Statement	Not applicable.

BPA Statement	Not applicable.
DEHP Statement	Not applicable.
RoHS	Not applicable.
REACH Statement	Not applicable.
Storage Conditions	Store at Room Temperature.
Shelf Life	2 years.

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
This document has been produced electronically and is valid without a signature.