

VWR® Cleanroom Clogs for critical environments

- 01. AUTOCLAVABLE
- 02. USE IN ISO 4 TO 6
(CLASS 10 TO CLASS 1000)
CLEANROOMS
- 03. ANTISTATIC



Comfort & durability

VWR® CLEANROOM CLOGS FOR CRITICAL CLEANROOM ENVIRONMENTS

Cleanroom clogs provide a smooth surface for day long comfort. A special conical design prevents spilled liquids from entering the clog.

- Antistatic
- Autoclavable
- Non-slip sole (Tested IAW ASTM F2913)
- Suitable for use in ISO 4 to 6 (Class 10 to Class 1000) cleanroom
- Heel strap
- Anatomic heel basin
- Toe grabber
- Color: White



They have been designed with foot comfort in mind. Smooth surface and aeration holes on both sides provide day extended comfort. Special conical design prevents spilled liquids from entering inside the clog.

These shoes are suitable for autoclave sterilisation. They can also be cleaned in 100 °C with domestic detergents, liquid disinfecting solutions, EO, or formaldehyde.

Ideal for use with VWR cleanroom disposable or launderable socks.

Size	Colour	Pk	Cat. No.
XS (EU: 35/36)	White	1 pair	114-8980
S (EU: 37/38)	White	1 pair	114-8981
M (EU: 39/40)	White	1 pair	114-8982
L (EU: 41/42)	White	1 pair	114-8983
XL (EU: 43/44)	White	1 pair	114-8984
XXL (EU: 45/46)	White	1 pair	114-8985

CONTACT

For more information contact your local sales specialist:
vwr.com/cleanroom-contact