

DuraCal pH buffer 10.01, 250 ml

Specification Sheet (Part/REF # 238321)

Patented buffer solutions provide pH stability of up to five years. High buffering capacity provides rapid and stable calibration.

Product Specifications

pH Accuracy	± 0.02 pH
Package	250 mL
Stability Months (after date of manufacturing)	60
Certificate	DAkkS
pH Value	10.01
Storage Conditions	5 - 35 °C / 41 - 95 °F

Specifications are subject to change without notice

<https://www.hamiltoncompany.com/process-analytics/sensors/238321>

Spec. Version B