

# VWR<sup>®</sup> Fixed Angle General Purpose Centrifuge

- 01. IDEAL FOR RESEARCH AND  
CLINICAL APPLICATIONS
- 02. 8 X 15 ML OR 6 X 50 ML  
CAPACITY
- 03. PULSE OPERATION FOR  
QUICK SPINS



VWR® FIXED ANGLE GENERAL PURPOSE CENTRIFUGE

The VWR® Fixed Angle General Purpose Centrifuge is ideal for research or clinical applications. This compact centrifuge can accept either an 8 x 15mL fixed angle rotor or a 6 x 50mL fixed angle rotor. There are also adapters available for 5, 7, and 10mL tubes. The digital display allows users to view either RPM or RCF as well as time. Variable speed can be adjusted between 300-4500 RPMs and the rotors are driven by a brushless DC motor drive. Additional features include a secure lid lock and over-speed detection for user safety.



Specifications

|                        |                                        |
|------------------------|----------------------------------------|
| Speed Range            | 300 to 4500 rpm, increments of 100 rpm |
| RCF Range              | 100 to 2490 x g, increments of 100     |
| Maximum Capacity       | 12 x 10 mL, 8 x 15 mL or 6 x 50 mL     |
| Timer                  | 30 sec to 99 min or continuous         |
| Dimensions (W x D x H) | 11 x 14.3 x 10.4 in.                   |
| Weight                 | 13.2 lbs                               |
| Electrical             | 100V - 240V, 50/60 Hz, 3A              |

| Description                                                       | Cat. No.  |
|-------------------------------------------------------------------|-----------|
| Centrifuge with 8 x 15 mL Fixed Angle Rotor                       | 76019-132 |
| Centrifuge without Rotor                                          | 76018-988 |
| Accessories                                                       |           |
| 8 x 15 mL Fixed Angle Rotor                                       | 76018-982 |
| 6 x 50 mL Fixed Angle Rotor                                       | 76018-976 |
| 3, 5, 7, and 10 mL Tube Adaptor for 15 & 50 mL Fixed Angle Rotor* | 76018-984 |
| 15 mL Tube Adapter for 50 mL Fixed Angle Rotor                    | 76018-978 |
| 3, 5, and 7 mL Adapter Stopper for 15 mL Fixed Angle Rotor        | 76018-986 |
| 3, 5, and 7 mL Adapter Stopper for 50 mL Fixed Angle Rotor        | 76018-980 |

\*when using this tube adapter in the 50 mL rotor, you must also purchase adapter 76018-978