

Certificate of Conformity

VWR Catalogue Number	See List
Description	VWR® Pipette Tips, Super Low retention, Natural
Country of Origin	Manufactured in China
Date of Issue (yyyy-mm-dd)	2022-04-04

Natural colour tips with Grey/Red Colour rack

Catalogue Number	Description	Sterility
612-6743	10 µl Super Low retention, Natural, DNase/RNase-free, Non-Pyrogenic, without filter	Non Sterile
612-6747	10 µl Super Low retention, Natural, DNase/RNase-free, Non-Pyrogenic, without filter	Sterile
612-6741	10 µl Super Low retention, Natural, DNase/RNase-free, Non-Pyrogenic, with filter	Sterile
612-6738	10 µl Long, Super Low retention, Natural, DNase/RNase-free, Non-Pyrogenic, without filter	Non Sterile
612-6748	10 µl Long, Super Low retention, Natural, DNase/RNase-free, Non-Pyrogenic, without filter	Sterile
612-6742	10 µl Long, Super Low retention, Natural, DNase/RNase-free, Non-Pyrogenic, with filter	Sterile

Natural colour tips with Grey/Yellow colour rack

Catalogue Number	Description	Sterility
612-6756	20 µl Super Low retention, Natural, DNase/RNase-free, Non-Pyrogenic, with filter	Non Sterile
612-6758	20 µl Super Low retention, Natural, DNase/RNase-free, Non-Pyrogenic, with filter	Sterile
612-6735	100 µl Super Low retention, Natural, DNase/RNase-free, Non-Pyrogenic, with filter	Non Sterile
612-6734	100 µl Super Low retention, Natural, DNase/RNase-free, Non-Pyrogenic, with filter	Sterile

Catalogue Number	Description	Sterility
612-6754	200 µl Super Low retention, Natural, DNase/RNase-free, Non-Pyrogenic, without filter	Non Sterile
612-6751	200 µl Super Low retention, Natural, DNase/RNase-free, Non-Pyrogenic, with filter	Non Sterile
612-6749	200 µl Super Low retention, Natural, DNase/RNase-free, Non-Pyrogenic, without filter	Sterile
612-6752	200 µl Super Low retention, Natural, DNase/RNase-free, Non-Pyrogenic, with filter	Sterile
612-6760	300 µl Super Low retention, Natural, DNase/RNase-free, Non-Pyrogenic, with filter	Non Sterile
612-6762	300 µl Super Low retention, Natural, DNase/RNase-free, Non-Pyrogenic, with filter	Sterile

Natural colour tips with Grey /Blue colour rack

Catalogue Number	Description	Sterility
612-6729	1000 µl Super Low retention, Natural, DNase/RNase-free, Non-Pyrogenic, without filter	Sterile
612-6769	1000 µl Super Low retention, Natural, DNase/RNase-free, Non-Pyrogenic, with filter	Sterile
612-6730	1000 µl Long, Super Low retention, Natural, DNase/RNase-free, Non-Pyrogenic, with filter	Sterile

Quality System Compliance

Products are manufactured under the **ISO 9001:2015 & ISO 13485:2016** standard. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements.

Trace Metal: All materials are certified as metal-free.

PCR Inhibitor: Products are tested via PCR amplification and gel electrophoresis analyses.

Product Specifications

This product is Class I category device as defined by the FDA in 21CFR Parts 862-892.

Material	Tips: Polypropylene (PP). Rack: Polypropylene (PP). Filter: Polyolefin (PO).
Sterilization	Product labeled as sterile is EB (Electron Beam) irradiated and dose released upon ISO11137 recommended practices in effect at the time of validation. Products labeled sterile meets a minimum requirement of 10 ⁻⁶ SAL with a specified dose range of 15-30 kGy. (which is not applicable for non-sterile products)
ATP Assay	This product has been tested and is ATP free.
DNase & RNase Free	This product is free of any detectable DNase/RNase contamination.
Nucleic Acid	This product has been tested and is free of any detectable Nucleic Acid.
BSE/TSE	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
Cytotoxicity	Testing is conducted to quality all material resins using ISO 10993 standards for cytotoxicity and have been shown to be non-toxic.
Non-Pyrogenic Statement	The acceptance level for product is 0.06 EU/mL (TAL Gel Clot Method). Test passed.
Latex Statement	This product is not made from natural rubber latex.
BPA Statement	Bisphenol are not used in the manufacture of the raw material and are not expected to be present.
DEHP Statement	Not applicable.
RoHS	No substances (Lead, Cadmium, Mercury, Hexavalent Chromium (Cr6+), Poly Brominated Biphenyls (PBB), Poly Brominated Diphenyl ethers (PBDE), Bis(2-Ethyhexyl) Phthalate (DEHP), Benzyl Butyl Phthalate (BBP), Dibutyl Phthalate (DBP), Diisobutyl Phthalate (DIBP)) are used in manufacturing the raw materials and final product. No routinely analyse is performed.

REACH Statement	Not applicable.
Storage Conditions	Store in Room temperature.
Shelf Life	The shelf-life for sterile product is 3 years after production. (which is not applicable for non-sterile products)

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
This document has been produced electronically and is valid without a signature