

PRODUCT SPECIFICATION

Product Code:	1025/28SC 2000ml Jacketed Beaker-PTFE screw thread connectors included.
Brand:	PYREX
Standards:	Borosilicate glass Type 3.3 meets the requirements of the following standards: ISO3585 and DIN 12217 United States and European Pharmacopoeia Type 1 glass
Graduations:	Printed with product Branding, capacity indicator and graduations (Colour: White enamel)
Dimensions:	External Diameter = 130mm Overall Height = 280mm Inner chamber capacity 2000ml Inner Diameter = 106mm Inner Height = 260mm
Pack Quantity:	Quantity 1 per pack

This is to certify that the goods detailed above have been inspected and tested in accordance with our BS EN ISO9001 and BS EN ISO13485 quality systems, and conform in all respects to the acceptable quality levels (AQL) detailed therein and any relevant specifications or drawings, unless recorded below as concessionary goods.

Irradiated goods are certified as sterile and a copy of their certification of sterility is available.

Signed: P Burton

Quality Technician

Date: 12 Jan 2021