

PRODUCT SPECIFICATION

REF 721500 CEDERROTH EYE WASH STATION



1 GENERAL PRODUCT DESCRIPTION

Description

An Eye Wash station containing two bottles of Cederroth Eye Wash and plastic and textile adhesive bandages, intended to be used in first aid situations. The products are protected against dust and moisture by transparent lids. Space for additional products is provided behind a lid with first aid instructions.

The Eye Wash bottles are instantly opened by twisting in either direction. The Adhesive bandages (plasters) are placed in dispensers for easy access and theft-proof storage. A special key is provided for refill replacement.

Medical devices in the station are CE-labelled.

Storage information

Dry storage in room temperature. For maximum protection, keep in original package until use.



2 DETAILED PRODUCT INFORMATION

Content per consumer unit

<u>Product</u>	<u>Qty</u>	<u>REF</u> <u>(refill)</u>	<u>Product</u>	<u>Qty</u>	<u>REF</u> <u>(refill)</u>
Eye wash station (empty)	1	-	First Aid Instruction (label)	1	-
Cederroth Eye Wash, 500 ml (2-pack)	1	725200	Key to the dispensers	1	51616620
Salvequick Plastic Plasters (45-pack) in dispenser	1	6036			
Salvequick Textile Plasters (40-pack) in dispenser	1	6444			

Dimensions (w x d x h):	(station)	558 x 286 x 117 mm
Weight:	(full station)	3184 g

Material

First aid station:	Green parts:	Polystyrene, ABS-plastic
	Transparent lid:	Polycarbonate
	Sealing profiles:	EPDM rubber (synthetic rubber)
First aid products:	Refer to Product specifications for each product	

Shelf life: 5 years from production date for each product in the station.

3 LOGISTIC INFORMATION

	CU	SKU	Pallet
GTIN:	7310807215006	7310807215006	7312027215002
Content:	1 station	1 CU	11 SKU/layer x 3 layers
Size (w x l x h mm):	294 x 572 x 123	294 x 572 x 123	800 x 1200 x 1032
Gross weight:	3520 g	3520 g	141 kg
Packing material:	Carton 313 g	-	-

