


VWR[®] Universal Low Retention Pipet Tip sustainability

- 01. **LIGHTEST RACK ON THE MARKET**
- 02. **BASES ARE REUSABLE AS A 96-WELL RESERVOIR**
- 03. **FULLY RECYCLABLE PACKAGING**





VWR® UNIVERSAL LOW RETENTION PIPET TIP SUSTAINABILITY

VWR Universal Low Retention Pipet Tips offer a variety of features to help you meet your sustainability goals.

Lightest rack in the world – By reducing the amount of plastic going in, it decreases the amount that needs to be recycled after use and decreases weight during transit, lowering the overall carbon footprint.

Racks are clear – This reduces the amount of colorant and makes racks easier to recycle at recycling facilities.

Tips have a thinner shrink wrap – This reduces the amount of plastic.

Reusable base – You can reuse the base of each rack as a 96-well reservoir.

Greenest reload on the market with CleanPak – 100% plastic with no contaminating cardboard. Cardboard is much more difficult to recycle than plastic as well, making this option a much greener choice. You also have an option of reusing the resealable plastic bags when you're finished with the tips.

Less plastic used in tips – Flexible proximal end of tip allows for less plastic use than our competition while providing a better fit and ergonomics.

30% less electricity than traditional systems – New generation of electric and hybrid molding machines increase efficiency.

Fewer defects and less waste – Six Sigma lean programs and premium raw materials maximize product quality and manufacturing productivity.

100% of packaging is fully recyclable – Additionally, inkjet is used on lids instead of stickers being placed on racks.



Volume	Length	Cat. No.	Pack of	Case of
Pre-Sterilized Tips				
10 µL	31.24 mm (1.23")	76322-512	960	4800
10 µL	45.84 mm (1.81")	76322-514	960	4800
200 µL	50.75 mm (2.00")	76322-516	960	4800
250 µL	59.54 mm (2.34")	76322-518	960	4800
300 µL	59.54 mm (2.34")	76322-520	960	4800
1000 µL	81.82 mm (3.22")	76322-522	960	4800
1250 µL	98.40 mm (3.87")	76322-138	960	4800

Packaging: 10 × 96 Tips/Pack

Volume	Length	Cat. No.	Pack of	Case of
Non-Sterile Tips				
10 µL	31.24 mm (1.23")	76322-140	960	4800
10 µL	45.84 mm (1.81")	76322-142	960	4800
200 µL	50.75 mm (2.00")	76322-144	960	4800
250 µL	59.54 mm (2.34")	76322-146	960	4800

Packaging: 10 × 96 Tips/Pack



VWR.COM

Prices, product, and/or services details are current when published and subject to change without notice. | Certain products or services may be limited by country, federal, state, provincial, or local regulations. | VWR, part of Avantor, makes no claims or warranties concerning sustainable/green products. Any claims concerning sustainable/green products are the sole claims of the manufacturer and not those of VWR International, LLC and/or Avantor, Inc. or affiliates. Offers valid in countries listed above, void where prohibited by law or company policy, while supplies last. | Trademarks are owned by Avantor, Inc. or its affiliates, unless otherwise noted. | Visit vwr.com to view our privacy policy, trademark owners, and additional disclaimers. © 2020 Avantor, Inc. All rights reserved.