



TECHNICAL DATA SHEET

HYGIENE MASK 3 PLY – Blue – Ear loops 113-0716

CHARACTERISTICS:

Hygiene mask with ear loops are made of non-woven and melt-blown polypropylene material, commonly used in a clean environment. The hygiene masks are manufactured and packed in a stringent clean environment.

PRODUCTS DESCRIPTION: Face Mask ear loop, glass free, lint-free, low particulating, no exposed metal

PHYSICAL PROPERTIES:

Material:	
Front Layer:	non-woven Polypropylene (blue 20 gram/m ²)
Back Layer:	non-woven Polypropylene short fiber (white 25 gram/m ²)
Center:	Melt-blown Polypropylene (white 20 gram/m ²)
Ear loop:	Latex-free elastic (round)
Size:	
Length:	175mm +/- 2mm
Height:	95mm +/- 2mm
Nosepiece:	115mm +/- 2mm
Ear loop:	175mm +/- 2mm
Colors	Blue outer, white inner, white ear loop
Bacterial Filtration Efficiency:	> 90 %
Particle Filtration Efficiency:	> 90 %
Particle size:	0.1Micron

PACKING:

50 pieces/box, 20boxes/carton

QUALITY ASSURANCE CRITERIA:

Hygiene masks are guaranteed to meet 4.0 AQL level II sampling plan per MIL-STD-105E for the following defects criteria:

- | | | | |
|----|------------------------------------|----|-----------------------|
| 1. | Ultrasonically welded not strong | 2. | Irregular PP material |
| 3. | Dirty PP material | 4. | Nosepiece exposed |
| 5. | Physical dimensions non-compliance | | |