



CERTIFICATE OF ANALYSIS



SIZE	10L (2.64gal)	1gal (3.8L)	115L (30gal)	210L (55gal)
DECON #	6001	6002	6003	6055
FISHER #	04-355-30	04-355-37	04-355-9	04-355-40

Lot #0D053 Manufacturing Date April 2020

Recommended Retest date is two years from the date of manufacture

TEST	SPECIFICATION	UNIT	TEST RESULT
pH	13 - 14	N/A	14
Density	1.3 - 1.4	g/ml	1.3
Appearance	Straw to amber colored clear liquid	N/A	Pass

Reviewed and Approved by:

Date: April 2020

A handwritten signature in blue ink, appearing to read 'Karen Wend', is written over a light blue horizontal line.

Karen Wend
Manager, Quality and Technical Services

BSE/TSE Statement: There are no materials of animal origin in the raw materials or in the process aides used to manufacture any Decon products.