


Axygen[®] Pipette Tips

High quality pipette tips for all your liquid handling needs

Experience the latest innovative solutions in liquid handling, from sustainable pipette tip racks, to optimized low retention technology.





Axygen® Pipette Tips

Overview	2
Axygen MaxRecovery+™ Technology	3
Axygen Universal Pipette Tips	4
Axygen LTS Compatible Tips	8
Axygen HybridRack® Pipette Tips	10
Axygen MultiRack Tip System	11
Related Products	13

Overview

Axygen® universal pipette tips are designed to fit a wide variety of single- and multi-channel pipettors. Flexible walls and a series of internal sealing rings ensure a secure fit with less force required to load and eject.

The Axygen product offering now also includes the Axygen MultiRack System, which is our updated reload and racking system. In addition to the universal fit tips packaged in the MultiRack System, we now also offer LTS Compatible tips in the MultiRack System, and if you are worried about your carbon footprint, the Axygen HybridRack® is a perfect option for you.

Automated assembly – Ensures product is free from contamination.

ISO 13485 and ISO 14001 – Certified quality management systems supply stringent quality processes to ensure the highest quality finished good.

In-house testing – In-house QA lab allows product to be validated for fit, form, and function prior to release.

Brand you can trust – A long track record of innovation and value as one of the Corning trusted brands.

Broad product offering – Corning offers many different styles of pipette tips, including bulk, racked, filtered, reloads, Tip Station, and our environmentally friendly HybridRack. With so many options, you are sure to find a tip that works for you.

Low retention technology – Pipette tips are available with the Axygen MaxRecovery+™ technology, our advanced technology to maximize sample recovery.




Packaging Solutions



Bulk




Reload Stack



Racked



Stacked
Rack



HybridRack

Corning products are included in Avantor's Environmentally Preferable Products (EPP) program, designed to provide greater transparency into product sustainability attributes. This allows you to make more informed and sustainable purchasing decisions.



Axygen[®] MaxRecovery+[™] Technology


Because every drop counts.

Axygen MaxRecovery+ pipette tips are our newest innovation utilizing the most advanced low retention technology yet. Developed through extensive research and backed by Corning's legacy of innovation and intellectual property protection, the MaxRecovery+ technology sets a new standard for low retention. With enhanced performance, improved solvent compatibility, and no intentionally added PFAS, the new surface is designed to maximize sample recovery while meeting regulatory and environmental standards. Whether handling aqueous or solvent-based samples, researchers can rely on MaxRecovery+ products for precision, accuracy, and superior performance. Upgrade your workflow with Axxygen MaxRecovery+ tips and experience the low retention performance and innovation you've come to expect from Corning.

Bringing Science to Life. Sustainably.

Corning is paving the way for cleaner air, cleaner water, and healthier generations to come. For Corning Life Sciences this means bringing science to life sustainably and moving forward in our three focus areas – climate initiatives, design for sustainability, and stewardship.

Go to www.corning.com/LifeSciESG to learn more.



Learn more about MaxRecovery+ Technology at www.corning.com/MaxRecovery.

Axygen® Universal Pipette Tips

Redesigned for the modern researcher.

For more than 20 years, Axxygen Universal pipette tips have set the standard for quality, accuracy, and precision in laboratories worldwide. Now redesigned to meet every researchers needs, the next generation of Axxygen Universal pipette tips deliver improved pipettor compatibility, durable hinged racks, and a complete offering from 10 µL to 1200 µL. Featuring tips with extended length, wide bore designs for viscous solutions, and our proprietary Axxygen MaxRecovery+™ technology, Axxygen Universal pipette tips are sure to meet and exceed your liquid handling expectations.

- ▶ RNase-/DNase-free
- ▶ Human gDNA-free
- ▶ Nonpyrogenic
- ▶ USP Class VI polypropylene

Axygen Universal Filter Pipette Tips

Racked, Filtered

	Standard Surface VWR Corning Cat. No.	MaxRecovery+™ VWR Corning Cat. No.	Volume (µL)	Pack Type	Sterile	Quantity
	78692-602 PTF-300-R-S	78692-600 PTF-300-L-R-S	0.5 - 10	Fine point	Yes	96/rack, 10 racks/pack, 5 packs/case
	78692-622 PTXLF-10-R-S	78692-620 PTXLF-10-L-R-S	0.5 - 10	Extended length, graduated, fine point	Yes	96/rack, 10 racks/pack, 5 packs/case
	78692-598 PTF-20-R-S	78692-596 PTF-20-L-R-S	1 - 20	Beveled	Yes	96/rack, 10 racks/pack, 5 packs/case
	78692-584 PTF-100-R-S	78692-582 PTF-100-L-R-S	1 - 100	Beveled	Yes	96/rack, 10 racks/pack, 5 packs/case
	78692-592 PTF-200-R-S	78692-590 PTF-200-L-R-S	1 - 200	Graduated, fine point	Yes	96/rack, 10 racks/pack, 5 packs/case
	- -	78692-594 PTF-205-WB-L-R-S	1 - 200	Wide bore	Yes	96/rack, 10 racks/pack, 5 packs/case
	78692-606 PTF-350-R-S	78692-604 PTF-350-L-R-S	1 - 300	Graduated, fine point	Yes	96/rack, 10 racks/pack, 5 packs/case
	78692-578 PTF-1000-R-S	78692-576 PTF-1000-L-R-S	100 - 1000	Beveled, graduated	Yes	96/rack, 10 racks/pack, 5 packs/case
	- -	78692-580 PTF-1005-WB-L-R-S	100 - 1000	Wide bore	Yes	96/rack, 10 racks/pack, 5 packs/case
	78692-586 PTF-1200-C-R-S	78692-588 PTF-1200-L-R-S	100 - 1200	Fine point	Yes	96/rack, 8 racks/pack, 4 packs/case

Reload, Filtered

	78692-452 PRFLF-300-C-S	78692-450 PRFLF-300-C-L-S	0.5 - 10	Fine point	Yes	96/deck, 5 decks/pack, 10 packs/case
	78692-466 PRFLT-XLF-10-C-S	78692-464 PRFLT-XLF-10-C-L-S	0.5 - 10	Extended length, graduated, fine point	Yes	96/deck, 5 decks/pack, 10 packs/case
	78692-448 PRFLF-200-C-S	78692-446 PRFLF-200-C-L-S	1 - 200	Graduated, fine point	Yes	96/deck, 5 decks/pack, 10 packs/case
	78692-444 PRFLF-1000-C-S	78692-442 PRFLF-1000-C-L-S	100 - 1000	Beveled, graduated	Yes	96/deck, 5 decks/pack, 10 packs/case

Axygen® Universal Pipette Tips

10 µL Fine Point



Standard Surface VWR Corning Cat. No.	MaxRecovery+™ VWR Corning Cat. No.	Volume (µL)	Pack Type	Sterile	Quantity
78692-552 PT-300	78692-554 PT-300-L	0.5 - 10	Bulk	No	1,000/pack, 10 packs/case
78692-560 PT-300-R	78692-556 PT-300-L-R	0.5 - 10	Racked	No	96/rack, 10 racks/pack, 5 packs/case
78692-562 PT-300-R-S	78692-558 PT-300-L-R-S	0.5 - 10	Racked	Yes	96/rack, 10 racks/pack, 5 packs/case
78692-438 PRFL-300-C	78692-440 PRFL-300-C-L	0.5 - 10	Reload	No	96 tips/deck, 10 decks/pack, 5 packs/case

10 µL Extended Length, Graduated, Fine Point



78692-608 PTXL-10	78692-610 PTXL-10-L	0.5 - 20	Bulk	No	1,000/pack, 10 packs/case
78692-616 PTXL-10-R	78692-612 PTXL-10-L-R	0.5 - 20	Racked	No	96/rack, 10 racks/pack, 5 packs/case
78692-618 PTXL-10-R-S	78692-614 PTXL-10-L-R-S	0.5 - 20	Racked	Yes	96/rack, 10 racks/pack, 5 packs/case
78692-460 PRFL-TXL-10-C	78692-462 PRFL-TXL-10-C-L	0.5 - 20	Reload	No	96/rack, 10 racks/pack, 5 packs/case

200 µL Clear



78692-510 PT-200-C	78692-512 PT-200-C-L	1 - 200	Bulk	No	1,000/pack, 10 packs/case
78692-522 PT-200-C-R	78692-514 PT-200-C-L-R	1 - 200	Racked	No	96/rack, 10 racks/pack, 5 packs/case
78692-524 PT-200-C-R-S	78692-516 PT-200-C-L-R-S	1 - 200	Racked	Yes	96/rack, 10 racks/pack, 5 packs/case
78692-526 PT-200-C-STK	78692-518 PT-200-C-L-STK	1 - 200	Stack rack	No	200/deck, 5 decks/ pack, 10 packs/case
78692-528 PT-200-C-STK-S	78692-520 PT-200-C-L-STK-S	1 - 200	Stack rack	Yes	200/deck, 5 decks/ pack, 10 packs/case
78692-434 PRFL-200-C	78692-436 PRFL-200-C-L	1 - 200	Reload	No	96 deck, 10 decks/pack, 5 packs/case

200 µL Beveled, Yellow



78692-530 PT-200-Y	78692-532 PT-200-Y-L	1 - 200	Bulk	No	1,000/pack, 10 packs/case
78692-538 PT-200-Y-R	78692-534 PT-200-Y-L-R	1 - 200	Racked	No	96/rack, 10 racks/pack, 5 packs/case
78692-540 PT-200-Y-R-S	78692-536 PT-200-Y-L-R-S	1 - 200	Racked	Yes	96/rack, 10 racks/pack, 5 packs/case
78692-542 PT-200-Y-STK	78692-542 PT-200-Y-STK	1 - 200	Stack rack	No	200/deck, 5 decks/ pack, 10 packs/case

200 µL Wide Bore



-	78692-546 PT-205-WB-C-L	1-200	Bulk	No	1,000/pack, 10 packs/case
-	78692-548 PT-205-WB-L-R	1-200	Racked	No	96/rack, 10 racks/pack, 5 packs/case
-	78692-550 PT-205-WB-L-R-S	1-200	Racked	Yes	96/rack, 10 racks/pack, 5 packs/case

Axygen® Universal Pipette Tips

300 µL Graduated, Fine Point



Standard Surface VWR Corning Cat. No.	MaxRecovery+™ VWR Corning Cat. No.	Volume (µL)	Pack Type	Sterile	Quantity
78692-564 PT-350-C	78692-566 PT-350-C-L	1 - 300	Bulk	No	1,000/pack, 10 packs/case
78692-572 PT-350-C-R	78692-568 PT-350-C-L-R	1 - 300	Racked	No	96/rack, 10 racks/pack, 5 packs/case
78692-574 PT-350-C-R-S	78692-570 PT-350-C-L-R-S	1 - 300	Racked	Yes	96/rack, 10 racks/pack, 5 packs/case

1000 µL Graduated, Beveled, Clear



78692-480 PT-1000-C	78692-482 PT-1000-C-L	100 - 1000	Bulk	No	1,000/pack, 10 packs/case
78692-488 PT-1000-C-R	78692-484 PT-1000-C-L-R	100 - 1000	Racked	No	96/rack, 10 racks/pack, 5 packs/case
78692-490 PT-1000-C-R-S	78692-486 PT-1000-C-L-R-S	100 - 1000	Racked	Yes	96/rack, 10 racks/pack, 5 packs/case
78692-430 PRFL-1000-C	78692-432 PRFL-1000-C-L	100 - 1000	Reload	No	96/deck, 10 decks/ pack, 5 packs/case

1000 µL Graduated, Beveled, Blue



78692-468 PT-1000-B	78692-470 PT-1000-B-L	100 - 1000	Bulk	No	1,000/pack, 10 packs/case
78692-476 PT-1000-B-R	78692-472 PT-1000-B-L-R	100 - 1000	Racked	No	96/rack, 10 racks/pack, 5 packs/case
78692-478 PT-1000-B-R-S	78692-474 PT-1000-B-L-R-S	100 - 1000	Racked	Yes	96/rack, 10 racks/pack, 5 packs/case

1000 µL Wide Bore



-	78692-492 PT-1005-WB-C-L	100 - 1000	Bulk	No	1,000/pack, 10 packs/case
-	78692-494 PT-1005-WB-C-L-R	100 - 1000	Racked	No	96/rack, 10 racks/pack, 5 packs/case
-	78692-496 PT-1005-WB-L-R-S	100 - 1000	Racked	Yes	96/rack, 10 racks/pack, 5 packs/case

1200 µL Fine point



78692-498 PT-1200-C	78692-500 PT-1200-C-L	1 - 1200	Bulk	No	1,000/pack, 10 packs/case
78692-506 PT-1200-R	78692-502 PT-1200-L-R	1 - 1200	Racked	No	96/rack, 8 racks/pack, 4 packs/case
78692-508 PT-1200-R-S	78692-504 PT-1200-L-R-S	1 - 1200	Racked	Yes	96/rack, 8 racks/pack, 4 packs/case

Empty Racks for Reloads

VWR Corning Cat. No.	Description	Sterile	Quantity
78692-458 PRFL-R-300	Empty Rack for 10 µL	No	10/pack, 5/case
78692-456 PRFL-R-200	Empty Rack for 200 µL	No	10/pack, 5/case
78692-454 PRFL-R-1000	Empty Rack for 1000 µL	No	10/pack, 5/case

Axygen® Macrovolume Pipette Tips

VWR Corning Cat. No.	Description	Pack Type	Sterile	Quantity
89005-638 T-5000-C	5 mL, clear	Bulk	No	250 tips/pack, 10 packs/case
89005-641 T-10ML-C	10 mL, clear	Bulk	No	100 tips/pack, 20 packs/case

Axygen® LTS Compatible Tips



Axygen LTS Compatible Tips are ergonomically designed to reduce discomfort and optimized to work with Rainin LiteTouch pipettors, all with the quality you have come to expect from Corning. By reducing the force required to seal and eject, Axxygen LTS Compatible tips can help you avoid using excess force to obtain an effective seal between the pipettor and pipette tip. Axxygen LTS Compatible tips are available in your favorite packaging formats including Axxygen HybridRack®, a sustainable packaging choice. All LTS tips are composed of USP Class VI polypropylene and manufactured in an ISO 13485 and ISO 14001 facility.

- ▶ USP Class VI polypropylene
- ▶ Sterile and nonsterile formats available
- ▶ Nonpyrogenic
- ▶ RNase-/DNase-free
- ▶ Human gDNA-free
- ▶ PCR inhibition-free

Axygen LTS Compatible Tips

20 µL Extended Length, Graduated, Fine Point



Standard Surface VWR Corning Cat. No.	MaxRecovery+™ VWR Corning Cat. No.	Pack Type	Sterile	Quantity
77395-566 LTF-20-R-S	78114-172 LTF-20-L-R-S	Racked, filtered	Yes	96/rack, 10 racks/pack, 5 packs/case
77395-598 LTFER-20-R-S	78114-188 LTFER-20-L-R-S	HybridRack, filtered	Yes	96/rack, 10 racks/pack, 5 packs/case
77395-572 LT-20-R-S	78114-156 LT-20-L-R-S	Racked	Yes	96/rack, 10 racks/pack, 5 packs/case
77395-578 LT-20-R	78114-176 LT-20-L-R	Racked	No	96/rack, 10 racks/pack, 5 packs/case
77395-584 LT-20-RL	78114-180 LT-20-L-RL	Reload stack	No	96/deck, 10 decks/pack, 4 packs/case
77395-590 LT-20-TS	78114-166 LT-20-L-TS	Tip station	No	192/layer, 5 layers/pack, 5 packs/case
77395-592 LT-20-TS-S	78114-174 LT-20-L-TS-S	Tip station	Yes	192/layer, 5 layers/pack, 5 packs/case
78018-168 LT-20-C	78114-164 LT-20-L-C	Bulk	No	1,000/pack, 10 packs/case

200 µL Graduated, Fine Point



77395-568 LTF-200-R-S	78114-152 LTF-200-L-R-S	Racked, filtered	Yes	96/rack, 10 racks/pack, 5 packs/case
77395-600 LTFER-200-R-S	78114-154 LTFER-200-L-R-S	HybridRack, filtered	Yes	96/rack, 10 racks/pack, 5 packs/case
77395-574 LT-200-R-S	78114-190 LT-200-L-R-S	Racked	Yes	96/rack, 10 racks/pack, 5 packs/case
77395-580 LT-200-R	78114-162 LT-200-L-R	Racked	No	96/rack, 10 racks/pack, 5 packs/case
77395-586 LT-200-RL	78114-160 LT-200-L-RL	Reload stack	No	96/deck, 10 decks/pack, 4 packs/case
77395-594 LT-200-TS	78114-184 LT-200-L-TS	Tip station	No	192/layer, 5 layers/pack, 5 packs/case
77395-596 LT-200-TS-S	78114-186 LT-200-L-TS-S	Tip station	Yes	192/layer, 5 layers/pack, 5 packs/case
78018-168 LT-200-C	78114-150 LT-200-L-C	Bulk	No	1,000/pack, 10 packs/case

Axygen® LTS Compatible Tips

1000 µL Graduated, Fine Point



Standard Surface VWR Corning Cat. No.	MaxRecovery+™ VWR Corning Cat. No.	Pack Type	Sterile	Quantity
77395-570 LTF-1000-R-S	78114-178 LTF-1000-L-R-S	Racked, filtered	Yes	96/rack, 8 racks/pack, 5 packs/case
77395-602 LTFER-1000-R-S	78114-168 LTFER-1000-L-R-S	HybridRack, filtered	Yes	96/rack, 10 racks/pack, 5 packs/case
77395-576 LT-1000-R-S	78114-148 LT-1000-L-R-S	Racked	Yes	96/rack, 8 racks/pack, 5 packs/case
77395-582 LT-1000-R	78114-158 LT-1000-L-R	Racked	No	96/rack, 8 racks/pack, 5 packs/case
77395-588 LT-1000-RL	78114-170 LT-1000-L-RL	Reload stack	No	96/deck, 8 decks/pack, 4 packs/case
78018-166 LT-1000-C	78114-182 LT-1000-L-C	Bulk	No	1,000/pack, 10 packs/case

Axygen® HybridRack® Pipette Tips



Born from a desire to reduce the amount of waste happening in laboratories, Axxygen HybridRack pipette tips are an example of our innovative spirit leading to a more sustainable future. By replacing the traditional plastic rack and lid normally found in racked pipette tip products with a paperboard rack and lid, we can reduce the amount of plastic per case by 70%, or approximately two pounds per case of 200 µL tips. The HybridRack is strong enough to support multi-channel pipetting, and as soon as the rack is empty, it can be recycled using any curbside recycling program.

All HybridRack SKUs have been third-party audited to the My Green Lab® ACT® Ecolabel Standard, the premier eco-nutrition label for laboratory products.

- USP Class VI polypropylene
- Sterile
- Nonpyrogenic
- RNase-/DNase-free
- Human DNA-free
- PCR inhibitor-free

Axygen HybridRack Pipette Tips

10 µL Extended Length, Graduated, Fine Point

Standard Surface VWR Corning Cat. No.	MaxRecovery+™ VWR Corning Cat. No.	Filtered	Sterile	Quantity
77269-858 ERF-10XT-R-S	77269-852 ERF-10XT-L-R-S	Yes	Yes	96/rack, 10 racks/pack, 5 packs/case
77269-846 ER-10XT-R-S	-	No	Yes	96/rack, 10 racks/pack, 5 packs/case

20 µL, 100 µL, and 200 µL Graduated, Fine Point

77269-866 ERF-20-R-S	77269-864 ERF-20-L-R-S	Yes	Yes	96/rack, 10 racks/pack, 5 packs/case
77269-868 ERF-100-R-S	77269-860 ERF-100-L-R-S	Yes	Yes	96/rack, 10 racks/pack, 5 packs/case
77269-876 ER-200-R-S	-	No	Yes	96/rack, 10 racks/pack, 5 packs/case

200 µL and 300 µL Graduated, Fine Point

77269-870 ERF-200NX-R-S	77269-862 ERF-200NX-L-R-S	Yes	Yes	96/rack, 10 racks/pack, 5 packs/case
77269-878 ER-300NX-R-S	-	No	Yes	96/rack, 10 racks/pack, 5 packs/case

300 µL Graduated, Extended Length, Fine Point

77269-872 ERF-300XT-R-S	77269-854 ERF-300XT-L-R-S	Yes	Yes	96/rack, 10 racks/pack, 5 packs/case
77269-848 ER-300XT-R-S	-	No	Yes	96/rack, 10 racks/pack, 5 packs/case

1250 µL Graduated, Extended Length, Fine Point

77269-874 ERF-1000XT-R-S	77269-856 ERF-1000XT-L-R-S	Yes	Yes	96/rack, 10 racks/pack, 5 packs/case
77269-850 ER-1000XT-R-S	-	No	Yes	96/rack, 10 racks/pack, 5 packs/case

To learn more about the Corning HybridRack tips, visit vwr.com/corning.

		US
Axygen® 10uL HybridRack™ Tip Salt Lake City, Utah, United States SKU ER-10XT-R-S		
Environmental Impact Scale Decreasing Environmental Impact		
Manufacturing		
Manufacturing Impact Reduction	1.0	
Renewable Energy Use	No	
Responsible Chemical Management	1.0	
Shipping Impact	1.0	
Product Content	5.0	
Packaging Content	1.0	
User Impact		
Energy Consumption (kWh/day)	N/A	
Water Consumption (gallons/day)	N/A	
Product Lifetime	6.0	
End of Life		
Packaging	1.0	
Product	6.6	
Environmental Impact Factor:		22.6
Label Valid Through:		June 2027
		act.mygreenlab.org



Axygen® MultiRack Tip System



Engineered to fit most pipettors on the market, the Axxygen MultiRack Tip System is a complete pipette tip system offering in various packaging styles to meet every lab's needs. The MultiRack Tip System offers standard and extended-length sizes in volumes from 0.1 µL to 1250 µL. Racked and reload packaging styles feature preloaded tips in easy-load, color-coded deck inserts. With a small benchtop footprint and reusable rack, this system minimizes environmental impact and maximizes lab space. Axxygen MultiRack Tip System products are manufactured in facilities that meet ISO standards.

- ▶ Racked options available
- ▶ USP Class VI polypropylene
- ▶ Sterile
- ▶ Nonpyrogenic
- ▶ RNase-/DNase-free
- ▶ Human DNA-free
- ▶ PCR inhibitor-free

Axygen MultiRack Universal Filter Pipette Tips



Standard Surface VWR Corning Cat. No.	MaxRecovery+™ VWR Corning Cat. No.	Volume	Description	Quantity
76320-724 MRF-10XT-R-S	76320-750 MRF-10XT-L-R-S	10 µL	Extended length, graduated, fine point	96/rack, 10 racks/pack, 5 packs/case
76320-726 MRF-20-R-S	76320-752 MRF-20-L-R-S	20 µL	Graduated, fine point	96/rack, 10 racks/pack, 5 packs/case
76320-728 MRF-100-R-S	76320-754 MRF-100-L-R-S	100 µL	Graduated, fine point	96/rack, 10 racks/pack, 5 packs/case
76320-730 MRF-200NX-R-S	76320-736 MRF-200NX-L-R-S	200 µL	Graduated, fine point	96/rack, 10 racks/pack, 5 packs/case
76320-732 MRF-200XT-R-S	76320-738 MRF-200XT-L-R-S	200 µL	Extended length, graduated, fine point	96/rack, 8 racks/pack, 5 packs/case
76320-734 MRF-300XT-R-S	76320-740 MRF-300XT-L-R-S	300 µL	Extended length, graduated, fine point	96/rack, 8 racks/pack, 5 packs/case
76320-748 MRF-1000XT-R-S	76320-742 MRF-1000XT-L-R-S	1250 µL	Extended length, graduated, fine point	96/rack, 8 racks/pack, 5 packs/case

Axygen MultiRack Universal Pipette Tips


10 µL Extended Length, Graduated, Fine Point



Standard Surface VWR Corning Cat. No.	MaxRecovery+ VWR Corning Cat. No.	Pack Type	Sterile	Quantity
76320-720 MR-10XT	-	Bulk	No	1,000/pack, 10 packs/case
76320-714 MR-10XT-RL	-	Reload	No	96/deck, 10 decks/pack, 4 packs/case
76320-676 MR-10XT-R	-	Racked	No	96/rack, 10 racks/pack, 5 packs/case
76320-686 MR-10XT-R-S	76320-696 MR-10XT-L-R-S	Racked	Yes	96/rack, 10 racks/pack, 5 packs/case
76320-706 MR-10XT-TS	-	Tip station	No	192/layer, 5 layers/pack, 5 packs/case
76320-708 MR-10XT-TS-S	-	Tip station	Yes	192/layer, 5 layers/pack, 5 packs/case

Axygen® MultiRack Universal Pipette Tips

200 µL, Graduated, Fine Point



Standard Surface VWR Corning Cat. No.	MaxRecovery+™ VWR Corning Cat. No.	Pack Type	Sterile	Quantity
76320-672 MR-200	- -	Bulk	No	1,000/pack, 10 packs/case
76320-716 MR-200-RL	- -	Reload	No	96/deck, 10 decks/pack, 4 packs/case
76320-678 MR-200-R	- -	Racked	No	96/rack, 10 racks/pack, 5 packs/case
76320-688 MR-200-R-S	76320-698 MR-200-L-R-S	Racked	Yes	96/rack, 10 racks/pack, 5 packs/case
76320-710 MR-200-TS	- -	Tip station	No	192/layer, 5 layers/pack, 5 packs/case
76320-712 MR-200-TS-S	- -	Tip station	Yes	192/layer, 5 layers/pack, 5 packs/case

300 µL, Graduated, Fine Point



Standard				
76320-718 MR-300-RL	- -	Reload	No	96/deck, 10 decks/pack, 4 packs/case
76320-680 MR-300-R	- -	Racked	No	96/rack, 10 racks/pack, 5 packs/case
76320-690 MR-300-R-S	76320-700 MR-300-L-R-S	Racked	Yes	96/rack, 10 racks/pack, 5 packs/case

300 µL Graduated, Extended Length, Fine Point



76320-682 MR-300XT-R	- -	Racked	No	96/rack, 8 racks/pack, 5 packs/case
76320-692 MR-300XT-R-S	76320-702 MR-300XT-L-R-S	Racked	Yes	96/rack, 8 racks/pack, 5 packs/case

1250 µL Graduated, Extended Length, Fine Point




76320-674 MR-1000XT	- -	Bulk	No	768/pack, 10 packs/case
76320-722 MR-1000XT-RL	- -	Reload	No	96/deck, 8 decks/pack, 4 packs/case
76320-684 MR-1000XT-R	- -	Racked	No	96/rack, 8 racks/pack, 5 packs/case
76320-694 MR-1000XT-R-S	- MR-1000XT-L-R-S	Racked	Yes	96/rack, 8 racks/pack, 5 packs/case

Axygen MultiRack Pipet Tip System, Empty Racks

VWR Corning Cat. No.	Description	Sterile	Quantity
76320-744 MR-R-A	Empty racks for 10 µL extended length tips, 200 µL, 200 µL NX, 300 µL tips	No	1/pack, 5/case
76320-746 MR-R-B	Empty racks for extended length 200 µL 300 µL, and 1000 µL tips	No	1/pack, 5/case

Related Products



Axygen® SnapLock Microcentrifuge Tubes

VWR Corning Cat. No.	Description	Quantity
10011-680 MCT-060-C	0.6 mL SnapLock microcentrifuge tubes, MaxyClear	1,000/pack, 10 packs/case
10011-682 MCT-060-C-S	0.6 mL SnapLock microcentrifuge tubes, MaxyClear, sterile	100/bag, 5 bags/pack, 10 packs/case
10011-700 MCT-150-C	1.5 mL SnapLock microcentrifuge tubes, MaxyClear	5,000/pack, 10 packs/case
10011-702 MCT-150-C-S	1.5 mL SnapLock microcentrifuge tubes, MaxyClear, sterile	50/bag, 5 bags/pack, 10 packs/case
10011-722 MCT-175-C	1.7 mL SnapLock microcentrifuge tubes, MaxyClear	500/pack, 10 packs/case
10011-724 MCT-175-C-S	1.7 mL SnapLock microcentrifuge tubes, MaxyClear, sterile	50/bag, 5 bags/pack, 10 packs/case
10011-742 MCT-200-C	2.0 mL SnapLock microcentrifuge tubes, MaxyClear	5,000/pack, 10 packs/case
10011-744 MCT-200-C-S	2.0 mL SnapLock microcentrifuge tubes, MaxyClear, sterile	50/bag, 5 bags/pack, 10 packs/case
89425-665 MCT-500-C	5.0 mL SnapLock microcentrifuge tubes, MaxyClear	250/pack, 5 packs/case
89425-665 MCT-500-C-S	5.0 mL SnapLock microcentrifuge tubes, MaxyClear, sterile	250/pack, 5 packs/case

NOTE: Other options are available, including full cases of different colors, assorted colors in the same case, and capless tubes.

For more related products, visit vwr.com/corning.

Warranty/Disclaimer: Unless otherwise specified, all products are for research use or general laboratory use only.* Not intended for use in diagnostic or therapeutic procedures. Not for use in humans. These products are not intended to mitigate the presence of microorganisms on surfaces or in the environment, where such organisms can be deleterious to humans or the environment. Corning Life Sciences makes no claims regarding the performance of these products for clinical or diagnostic applications. *For a listing of US medical devices, regulatory classifications or specific information on claims, visit www.corning.com/resources.

Corning's products are not specifically designed and tested for diagnostic testing. Many Corning products, though not specific for diagnostic testing, can be used in the workflow and preparation of the test at the customers discretion. Customers may use these products to support their claims. We cannot make any claims or statements that our products are approved for diagnostic testing either directly or indirectly. The customer is responsible for any testing, validation, and/or regulatory submissions that may be required to support the safety and efficacy of their intended application.

CORNING

For additional product or technical information, visit vwr.com/corning, call 1.800.932.5000, or contact your VWR representative.



vwr.com

Any prices, product, and/or services details are current when published and subject to change. Offers subject to product availability. Void where prohibited by law or company policy. Any claims concerning products manufactured by third parties are the sole claims of the manufacturer and not those of Avantor or its affiliates. All products are subject to the terms and conditions of sale of Avantor or its affiliates, including the applicable limitations of liability contained therein. The information herein should not be considered a warranty of any kind. Please carefully review any product information sheets and safety information for the product. Trademarks are owned by Avantor Inc. or its affiliates, unless otherwise noted. | Visit vwr.com to view our privacy policy and additional disclaimers. © 2026 Avantor Inc. All rights reserved.