

Certificate of Conformity

VWR Catalogue Number	113-9756
Description	Cleanroom Face Mask with ear loops, Class 100, Colour White, Length 17.5 ± 0.5 cm, Width 9.5 ± 0.2 cm
Lot Number	201116
Date of Manufacture (yyyy-mm-dd)	2020-11-16
Expiry Date (yyyy-mm-dd)	2023-11-16
Country of Origin	Manufactured in China
Date of Issue (yyyy-mm-dd)	2021-08-16

Quality System Compliance

Products are manufactured under the **ISO 9001:2015 & ISO 13485:2016** standard. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable product specifications/industry standards EN 14683:2019 TYPE IIR. Products are guaranteed to meet 4.0 AQL level II.

	Test / Standard	Specification	Result
Bacterial Filtration Efficiency [BFE]	ASTM F2101	> 95%	Pass
Pressure Drop [Delta P] [mmH ₂ O]	ASTM F2101	< 6	Pass

Product Specifications

Material	18 gsm ES + 25 gsm MB + 18 gsm ES.
Sterilization	Not Applicable.
ATP Assay	Not Applicable.
DNase & RNase Free	Not Applicable.

BSE/TSE	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
Cytotoxicity	Not Applicable.
Non-Pyrogenic Statement	Not Applicable.
Latex Statement	This product is not made from natural rubber latex.
BPA Statement	Not Applicable.
DEHP Statement	Not Applicable.
RoHS	Not Applicable.
REACH Statement	Not Applicable.
Storage Conditions	Room temperature, dry and shady.
Shelf Life	3 years.

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
This document has been produced electronically and is valid without a signature.