

Technical Datasheet

VWR WIPES SUPERSOFT

VWR International certifies that the products comply in every respect with the relevant particular specifications and/or drawings and/or relevant standards and regulations in force, mentioned in this document.

VWR International certifie que les produits sont conformes à tous égards aux spécifications particulières appropriées et / ou dessins, normes et règlements en vigueur, mentionnées dans le présent document

VWR International bestätigt, dass das Produkte allen zutreffenden Spezifikations- und/oder Zeichnungs- und/oder Standardanforderungen sowie den im Dokument aufgeführten gültigen Vorschriften entspricht.

Applicable part numbers : 115-2416, 115-2417

Attribute	Value*	Test method
Composition :		
viscose	60%	
polyester	40%	
Production method	spunlace hydrojet	
Colour	white	
Basis weight	50 g/m ²	ERT 40.3-90
Absorption	1010%	ERT 10.4-02
Tensile strength :	Dry Wet	ERT 20.2-89
Machine direction	94,0 N 74,7 N	
Cross direction	22,1 N 19,2 N	

Characteristics :

Extremely soft ; Non-scratching ; Very absorbent ; Resistant to bases and solvents ; Hygienic.
Recyclable. Does not contain binders.

Applications :

Automotive, High-Tech, Optics, Industry, Laboratories, Healthcare, Janitorial.

Formats & Packaging :

Cat No	Type	Format size	Quantity	Primary pkg	Secondary pkg	Pallet
115-2416	Industrial roll	40x30cm	2 rolls x 375 pcs	Plastic film	n/a	24 x 2rolls
115-2417	Z-folded wipes	38x32cm	10 packs x 50 pcs	Plastic bag	Cardboard case	48 cases

*** Disclaimer:** The information presented here is the result of testing and other evaluations conducted during the development of the product. This information does not constitute a specification.

Date of issue: **March 2017**