

Specification

Sodium Azide pure

A1430

Physical Description:	Solid
Product Code:	A1430
Product Name:	Sodium Azide pure
Specifications:	Assay (titr.): min. 99 % H₂O-insoluble matter: max. 0.05 % Solubility (H₂O; 17°C): 420 g/L Loss on drying: max. 0.25 % Carbonate (as Na₂CO₃): max. 0.25 %
Hazard pictograms	
UN:	1687
Class/PG:	6.1/II
ADR:	6.1/II
IMDG:	6.1/II
IATA:	6.1/II
WGK:	2
Storage:	RT
Signal Word:	Danger
GHS Symbols:	GHS06 GHS09
H Phrases:	EUH032 H300+H310 H410

P Phrases:

P280

P301+P310

AppliChem GmbH

Ottoweg 4 • D-64291 Darmstadt • Phone +49 6151 9357 0 • Fax +49 6151 9357 11 • info.de@itwreagents.com • www.itwreagents.com
 CEO Joan Roget • Commerzbank Darmstadt • Bank 508 800 50 • Account 0186989900 IBAN DE24 5088 0050 0186 9899 00 • Swiftcode DRESDEFF508 • Finanzamt Darmstadt 0072250192 • Register court Darmstadt HRB Nr. 7340

Specification

Sodium Azide pure

A1430

	P321
	P330
	P361+P364
	P405
	P501
Molecular Formula:	NaN ₃
M:	65.01 g/mol
CAS:	26628-22-8
EINECS:	247-852-1
CS:	28500060
Index Nr.:	011-004-00-7

AppliChem GmbH

Ottoweg 4 • D-64291 Darmstadt • Phone +49 6151 9357 0 • Fax +49 6151 9357 11 • info.de@itwreagents.com • www.itwreagents.com
CEO Joan Roget • Commerzbank Darmstadt • Bank 508 800 50 • Account 0186989900 IBAN DE24 5088 0050 0186 9899 00 • Swiftcode DRESDEFF508 • Finanzamt Darmstadt 0072250192 • Register court Darmstadt HRB Nr. 7340