

Certificate of Conformity

VWR Catalogue Number	VWRA548-2655
Description	VWR® Inserts for Screw Neck Vials ND10, Micro-insert with mandrel interior and polymer feet, suits for 10-425 vials, 250 µL volume, 29 × 6 mm dimension
Lot Number	20250410
Date of Manufacture (yyyy-mm-dd)	2025-04-10
Country of Origin	Manufactured in China
Date of Issue (yyyy-mm-dd)	2025-08-12

Quality System Compliance

This certificate confirms that the product has been sourced from a supplier or manufacturer who has declared their Quality Management System (QMS) complies with [GB/T 19001-2016/ISO 9001:2015](#).

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements.

Test/Procedure	Result
Inside Surface Hydrolytic Resistance	Accord with HC2
Coefficient of Thermal Expansion	(7.6~9.0)*10 ⁻⁶ K ⁻¹ (20°C~300°C), Accord with level 2.
Hydrolytic resistance of glass grains at 121 °C	Accord with level 2
Thermal shock endurance resistance	Pass
Stress	Pass
Leaching amount of As、Sb、Pb、Cd	As≤0.2mg/L ; Sb≤0.7mg/L ; Pb≤1.0mg/L ; Cd≤0.25mg/L.

Product Specifications

Material	(7.6~9.0) Low Borosilicate glass, Colour: Clear
Sterilization	Not applicable.
ATP Assay	Not applicable.

DNase & RNase Free	Not applicable.
BSE/TSE	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE) according to EMEA 410/01 revision 3.
Cytotoxicity	Not applicable.
Non-Pyrogenic Statement	Not applicable.
Latex Statement	The product is Latex free.
BPA Statement	Not applicable.
DEHP Statement	Not applicable.
RoHS	Raw materials comply with the requirements of the Directive 2002/95/EC, as amended, and 2011/65/Eu concerning the limits of cadmium, lead, mercury, hexavalent chromium, Polybrominated Biphenyls (PBB) and Polybrominated Diphenyl Ethers (PBDE).
REACH Statement	Not applicable.
Storage Conditions	Room temperature. away from direct sunlight and rain.
Shelf Life	Not applicable.

Disclaimer: VWR's Quality Management System (QMS) is ISO 9001:2015 certified as a wholesale distributor of finished products, ensuring effective oversight of finished product handling, traceability, and distribution processes. VWR is not a manufacturer and does not engage in manufacturing, co-manufacturing or processing activities such as design, packaging, labeling, testing or quality release. Consequently, the supplier or manufacturer's QMS certification may differ from VWR's.

VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. This product is intended for research use only. It is not for use in diagnostic procedures, therapeutic applications, or any human or clinical use.

This document has been produced electronically and is valid without a signature.