



Certificate of Conformity

VWR WINGED PLUG BL LDPE 13mm

VWR International certifies that the specified products meet all VWR specifications and quality criteria for release. The products comply in every respect with the relevant particular specifications, drawings, relevant standards and regulations in force.

VWR International certifie que les produits mentionnés respectent toutes les spécifications et critères de qualité VWR. Les produits satisfont aux spécifications particulières, aux plans, aux normes et règlements en vigueur s'y rapportant.

VWR International bestätigt hiermit, daß die angegebenen Produkte alle VWR Spezifikationen und Qualitätskriterien zur Freigabe erfüllen. Die Produkte entsprechen in jeder Hinsicht den maßgeblichen Spezifikationen, den Zeichnungen sowie den relevanten Normen und gültigen Vorschriften.

VWR Catalog No: **212-0577**

Cap Material: **Polyethylene**

Appearance: **Cylindrical**

Outside sizes: **ø : 13 mm**

Height cap: **18 mm**

Colour: Cap: **Blue**

Box Quantity: **20.000 pcs box**

Box Volume: **52.5 x 37.5 x 35 cm**

Box Weight: **13.5 kg**

Storage: **Do not expose or leave near open flames, do not leave the product around as it is not biodegradable.**

Disclaimer: VWR wants to state that these specifications can only be regarded as general information and will not discharge the user from his obligation to ensure the product is suitable for the intended use - If there is any doubt ask us.

Date of issue: **March 2015**