



The square COMBI BOX is equipped with a unique safety closure system that guarantees total leak-proofness. It effectively protects sensitive samples, solutions, or preparations from leakage or contamination. The tamper-proof lid prevents accidental opening and ensures the integrity of the contents until they are used, an essential advantage in research, analysis, or liquid storage environments.

Made from food-grade polypropylene (PP), COMBI BOX boxes offer perfect transparency for immediate visual inspection of stored samples and solutions. This visibility facilitates rapid identification of contents without having to open the boxes, thereby reducing unnecessary handling. Their square design optimizes organization in cabinets, refrigerators, or storage areas, while allowing for stable and secure stacking.

PRODUCT ADVANTAGES

- **Total leak-proof** : prevents any liquid leakage, even for fluid solutions.
- **Tamper-proof closure** : secures samples and guarantees the integrity of the contents.
- **Contamination prevention** : limits the risks associated with accidental opening.
- **Optimal transparency** : immediate visual inspection of contents without opening.
- **High-quality polypropylene (PP)** : inert material, resistant to common chemicals and sterilizable.

SPECIFICATIONS

SKU	MULC47534
Ext. dim.	193,28 x 193,28 x 76,95 mm
Volume	2L
Material	Polypropylene
Type of walls / base	Solid / solid
Color	Transparent
Food Safe	Yes
Stackable / nestable	Yes
Recyclable	Yes
Lid	Sold with lid
Shipping time	5 working days