

Certificate of Conformity

VWR Catalogue Number	231-0593
Description	Disposable Spatulas, Sterile, Non pyrogenic, 210 mm length, 7 mm blade width, Natural colour
Lot Number	031323
Date of Manufacture (yyyy-mm)	2023-03
Date of Sterilization (yyyy-mm)	2023-05
Expiry Date (yyyy-mm)	72 MONTHS
Country of Origin	Manufactured in USA
Date of Issue (yyyy-mm-dd)	2023-08-17

Quality System Compliance

Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements.

Product Specifications

Material	100 % Food Grade Polypropylene Homopolymer. Product complies with the FDA Regulations 21CFR 177.1520.
Sterilization	The products is sterilized to an SAL 10 ⁻⁶ using ethylene oxide as an Sterilant in compliance with ISO 11135 and ISO 11607-1 requirements. Autoclave: Autoclaving polypropylene at 121°C at 15 psi for 20minutes is standard.
ATP Assay	Not Applicable.
DNase & RNase Free	Not Applicable.

BSE/TSE	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
Cytotoxicity	Testing is conducted to quality all material resins using ISO 10993-7 standards for cytotoxicity and found to be EO <4 mg/device, ECH <9 mg/device.
Non-Pyrogenic Statement	The acceptance level for the non pyrogenic product is ≤ 0.5 EU/ml or 20 EU/device per USP <85> requirements.
Latex Statement	This product is not made from natural rubber latex.
BPA Statement	Bisphenol are not used in the manufacture of the raw material and are not expected to be present.
DEHP Statement	The product does not contain Di(2-Ethyhexyl) Phthalate (DEHP).
RoHS	Not Applicable.
REACH Statement	Not Applicable.
Storage Conditions	Storage of the spatulas out of direct sunlight in a cool, dry locations (normal lab conditions).

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
This document has been produced electronically and is valid without a signature.