

Certificate of Conformity

VWR Catalogue Number	See List
Description	VWR®, Cryo Boxes
Country of Origin	Manufactured in India
Date of Issue (yyyy-mm-dd)	2023-12-29

Catalogue Number	Compartments	Colour	For Tubes [mL]	W x D x H [mm]
479-2388	25	Assorted	1 - 2	75.6 × 75.6 × 63.0
479-2389	25	Blue	1 - 2	75.6 × 75.6 × 63.0
479-2390	25	Green	1 - 2	75.6 × 75.6 × 63.0
479-2391	25	Red	1 - 2	75.6 × 75.6 × 63.0
479-2392	25	White	1 - 2	75.6 × 75.6 × 63.0
479-2393	25	Yellow	1 - 2	75.6 × 75.6 × 63.0
479-2394	81	Assorted	1 - 2	132× 132 × 158
479-2395	81	Blue	1 - 2	132× 132 × 158
479-2396	81	Green	1 - 2	132× 132 × 158
479-2397	81	Red	1 - 2	132× 132 × 158
479-2398	81	White	1 - 2	132× 132 × 158
479-2399	81	Yellow	1 - 2	132× 132 × 158
479-2400	81	Assorted	3	132 × 132 × 179.5
479-2401	81	Blue	3	132 × 132 × 179.5
479-2402	81	Green	3	132 × 132 × 179.5
479-2403	81	Red	3	132 × 132 × 179.5
479-2404	81	White	3	132 × 132 × 179.5
479-2405	81	Yellow	3	132 × 132 × 179.5
479-2412	81	Assorted	4 - 5	132 × 132 × 200
479-2414	81	Blue	4 - 5	132 × 132 × 200
479-2415	81	Green	4 - 5	132 × 132 × 200

Catalogue Number	Compartments	Colour	For Tubes [mL]	W x D x H [mm]
479-2416	81	Red	4 - 5	132 x 132 x 200
479-2417	81	White	4 - 5	132 x 132 x 200
479-2418	81	Yellow	4 - 5	132 x 132 x 200
479-2406	100	Assorted	1 - 2	132 x 132 x 170
479-2407	100	Blue	1 - 2	132 x 132 x 170
479-2408	100	Green	1 - 2	132 x 132 x 170
479-2409	100	Red	1 - 2	132 x 132 x 170
479-2410	100	White	1 - 2	132 x 132 x 170
479-2411	100	Yellow	1 - 2	132 x 132 x 170

Quality System Compliance

Products are manufactured under the **ISO 9001:2015 & ISO 13485:2016** standard. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements.

Test Name	Test Reference	Accuracy	Result
Autoclaving	QWI-QC-07	General Insp. Level-1, AQL -1.5	PASS
Dimensional	QWI-QC-23	General Insp. Level-1, AQL -1.5	Within tolerance

Product Specifications


Material	Polycarbonate.
Sterilization	Non-sterile.
ATP Assay	Not Applicable.
DNase & RNase Free	Not Applicable.
BSE/TSE	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).

Cytotoxicity	Not Applicable.
Non-Pyrogenic Statement	Not Applicable.
Latex Statement	The products are not made from natural rubber latex.
BPA Statement	Not Applicable.
DEHP Statement	Not Applicable.
RoHS	Not Applicable.
REACH Statement	Not Applicable.
Storage Conditions	Application temperature -196 °C to +121 °C. Product to be kept in a cool, dry rack at room temperature. The product should not be exposed to the sunlight or extreme temperatures.
Shelf Life	Not Applicable.

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
This document has been produced electronically and is valid without a signature.

Annex

Technical Drawing



Version History

Sl.no	Version number	Description	Month & Year
1	1.0	Created General CoC using new layout	12-2023