


	TECHNICAL DATA SHEET	400400	Date: 11/05/16
		Revision 1	Page 1/3

**Codes: 400500, 400800, 400800.1, 400700, 400400, 401700,
400900, 400705, 401100**
Ultra clear polypropylene tubes

General Catalogue page 122

Mod.	Code	Dimensions (mm)	Volume (ml)	Case Quantity	Case weight (kg)	Case volume	Cases per pallet
1	400500	11 x 55	3	2 x 2,000	6.50	0.046	32
2	400800	12 x 75	5	4 x 1,000	8.00	0.075	20
3	400800.1	12 x 75	5	4 x 1,000	8.00	0.075	20
4	400700	13 x 75	5	4 x 1,000	9.70	0.082	20
5	400400	12 x 88	6	4 x 1,000	9.60	0.082	20
6	401700	13 x 100	7	6 x 500	9.40	0.082	20
7	400900	16 x 95	10	5 x 500	10.10	0.096	16
8	400705	16 x 100	10	5 x 500	12.20	0.096	16
9	401100	15 x 50	5	5 x 1,000	10.40	0.082	20

Directive 98/79/EC. "In vitro" diagnostic medical devices



Written	Verified	Authorised
 Gemma Riba Marketing Dept.	 Anna Mir Technical, Development and Quality Director	 Monica Torras Sales & Marketing Dept.

	TECHNICAL DATA SHEET	400400	Date: 11/05/16
		Revision 1	Page 2/3

USE AND DESCRIPTION

Manufactured from high technology, new generation moulds to ensure the reproducibility of each tube. Tubes are made in one piece, without joint, to ensure uniformity and dimensional accuracy from tube to tube. All tubes are suitable for laboratory use.

All tubes are suitable for centrifugation up to **15,000 xg**

The following codes have graduation rings:

- **400800** (mod. 2): graduated at 1, 2.5 and 4 ml
- **400400** (mod. 5): graduated at 1, 2.5 and 5 ml
- **400900** (mod. 7): graduated at 2.5, 5 and 10 ml
- **401100** (mod. 9): graduated at 2.5 and 5 ml

MEASURES

Code	Internal Ø (mm)	External Ø (mm)	Total length (mm)	Maximum volume (ml)
400500	9.9	11.4	53.9	3.4
400800	10.3	11.7	73.9	5.3
400800.1	10.2	11.7	73.9	5.3
400700	10.9	12.7	73.7	6.4
400400	10.3	11.9	86.9	5.8
401700	11.0	12.7	98.4	8.1
400900	13.7	15.5	93.7	11.9
400705	13.7	15.8	98.6	12.7
401100	13.7	15.5	48.9	6.5

Written	Verified	Authorised
 Gemma Riba Marketing Dept.	 Anna Mir Technical, Development and Quality Director	 Monica Torras Sales & Marketing Dept.

	TECHNICAL DATA SHEET	400400	Date: 11/05/16
		Revision 1	Page 3/3

SUITABLE CAPS (NOT INCLUDED)

Suitable caps:

- **400500**: code 3053
- **400800**: code 3058
- **400800.1**: code 3058
- **400700**: code 3071
- **400400**: code 3058
- **401700**: code 3071
- **400900**: codes 3072, 3050 and 3081
- **400705**: codes 3072, 3050 and 3081
- **401100**: codes 3072, 3050 and 3081

Colours (add number at the end of the code):

- **Natural (02)**: codes 3053, 3058, 3071, 3072, 3050, 3081
- **Yellow (01)**: code 3071
- **Malve (08)**: codes 3053, 3058, 3071, 3058
- **Blue (05)**: codes 3058, 3071, 3072
- **Green (07)**: codes 3053, 3058, 3071, 3072, 3050, 3081
- **Red (06)**: codes 3053, 3071, 3072, 3050, 3081
- **Black (09)**: codes 3058, 3071



Written	Verified	Authorised
 Gemma Riba Marketing Dept.	 Anna Mir Technical, Development and Quality Director	 Monica Torras Sales & Marketing Dept.