

# Certificate of Conformity

|                                   |                                          |
|-----------------------------------|------------------------------------------|
| <b>VWR Catalogue Number</b>       | See List                                 |
| <b>Description</b>                | Biohazard specimen bags, 50 µm thickness |
| <b>Country of Origin</b>          | Manufactured in USA                      |
| <b>Date of Issue (yyyy-mm-dd)</b> | 2024-03-29                               |

| Catalogue Number | Material | Length x Width [mm] |
|------------------|----------|---------------------|
| 129-1576         | LDPE     | 228 × 152           |
| 129-1577         | LDPE     | 304 × 228           |

## Quality System Compliance

Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

## QC Testing

Representative products samples are collected and inspected in accordance with internationally recognized ASTM Test methods for Dart Impact Strength (ASTM D1709) and tear resistance (ASTM D1922).

## Product Specifications

|                               |                                |
|-------------------------------|--------------------------------|
| <b>Material</b>               | Polyethylene Polymer (PE).     |
| <b>Sterilization</b>          | Non-Sterile, Non-autoclavable. |
| <b>ATP Assay</b>              | Not Applicable.                |
| <b>DNase &amp; RNase Free</b> | Not Applicable.                |

|                                |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
|--------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>BSE/TSE</b>                 | The film and resin suppliers confirm that some films, resins and masterbatch additive packages used to manufacture polyethylene films, may contain byproducts of bovine sourced raw materials in their manufacture. As advised by these suppliers, these animal byproduct additives are processed under highly refined and controlled industrial chemical processes and at very high temperatures prior to purification and subsequent use as a polymer additive. This applies to all products comprised of either polypropylene or low-density polyethylene polymers. The animal derived byproducts are not added during the process of manufacture of the product but originate either in the manufacture of the resins or in the additives used in film formulations. |
| <b>Cytotoxicity</b>            | Not Applicable.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| <b>Non-Pyrogenic Statement</b> | Not Applicable.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| <b>Latex Statement</b>         | This product is not made from natural rubber latex.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| <b>BPA Statement</b>           | Not Applicable.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| <b>DEHP Statement</b>          | Please see RoHS information below.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| <b>RoHS</b>                    | No substances (Lead, Cadmium, Mercury, Hexavalent Chromium (Cr6+), Poly Brominated Biphenyls (PBB), Poly Brominated Diphenyl ethers (PBDE), Bis(2-Ethyhexyl) Phthalate (DEHP), Benzyl Butyl Phthalate (BBP), Dibutyl Phthalate (DBP), Diisobutyl Phthalate (DIBP)) are used in manufacturing the raw materials and final product. No routinely analyse is performed.                                                                                                                                                                                                                                                                                                                                                                                                     |
| <b>REACH Statement</b>         | No substances as listed on >Substance of Very High Concern< candidate list published 2024-01-23 are used in manufacturing of the raw materials used to produce this product. No routinely analysis is performed.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| <b>Storage Conditions</b>      | Product may react with strong oxidizing agents and should not be stored near such materials. Store boxes and bags of material away from sunlight and in dry areas protected with automatic sprinklers. Storage temperature: 60° C (140° F) maximum.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| <b>Shelf Life</b>              | Not applicable.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.  
This document has been produced electronically and is valid without a signature.

## Version History

| Sl.no | Version number | Description                          | Month & Year |
|-------|----------------|--------------------------------------|--------------|
| 1     | 1.0            | Created General CoQ using new layout | 03-2024      |
|       |                |                                      |              |
|       |                |                                      |              |